

# TOMMANIKIN SIMULATOR RACK

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



The Tommanikin® Simulator Rack is a 3-shelf storage rack made of welded steel. It is designed to accommodate three Tommanikin® on NAR® Litters (Talon®, Raven®, T3 and MedEvac), as well as other NATO Litters (sold separately). Four locking casters are included with this product.

*Please note that this product is not intended to be used for patient transport*

## FEATURES/SPECIFICATIONS

- Durable, made of welded steel
- Can store up to three mannequins and litters (sold separately)
- Four locking casters included



**Litters and TOMMANIKIN® simulators not included**



## TOMMANIKIN RACK

ITEM #	NSN
93-0135	----

*\* Explore more Techline Trauma products at [www.NARescue.com](http://www.NARescue.com).*

## DIMENSIONS

- Dimensions: H 36 in. x W 36 in. x D 26 in.
- Weight: 150 lbs.