

BLEEDING CONTROL SKILLS TRAINING KITS

Learn more about this and other Products with a Mission® at www.NARescue.com

WITH AND WITHOUT WOUND MANAGEMENT SIMULATOR



BE PREPARED to save lives & stop life-threatening bleeding in the first critical minutes!

The Bleeding Control Skills Training Kits from North American Rescue® provide the appropriate equipment needed to teach up to 20 students the necessary life-saving skills of bleeding control. The equipment is packaged in a self-contained, wheeled rigid case with internal dividers to keep the equipment organized and easy to inventory. The instructor kit is available with or without a wound management simulator model. The contents of this advanced training kit include all the necessary equipment to address extremity bleeding, wound packing, pressure bandages, as well as litters to carry injured victims, corresponding to the equipment in NAR's Basic / Advanced / Advanced BCD / Advanced CELOX™ Public Access Bleeding Control Kits.



Available with or without Wound Management Simulator

FEATURES/SPECIFICATIONS

- Provides the appropriate equipment to teach up to 20 students the necessary life-saving skills of bleeding control
- Packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory
- Skill set training addresses extremity bleeding, wound packing, pressure bandages, chest injuries and patient transport
- Corresponds with NAR's Basic / Advanced / Advanced BCD / Advanced Celox™ Public Access Bleeding Control Kits

HARD CASE SPECIFICATIONS

- Retractable extension handle
- 4 strong polyurethane wheels with stainless steel bearings
- Easy-open double throw latches
- Open cell core with solid wall design – strong & lightweight
- O-ring seal
- Automatic pressure equalization valve
- Fold down handles
- Stainless steel hardware and padlock protectors
- Black color

DIMENSIONS

- Hard Storage Case: W 32 in.x H 21 in. x D 12.5 in.
- Total Weight
 - » With Wound Mgmt Sim: 48 lb
 - » Without Wound Mgmt Sim: 42 lb

BLEEDING CONTROL SKILLS TRAINING KITS

ITEM #	OPTION	WOUND MGMT SIMULATOR
80-0880	Basic	with WMS
80-0878	Intermediate	with WMS
80-0876	Advanced	with WMS
80-1584	Advanced with CELOX™	with WMS
80-0881	Basic	no WMS
80-0879	Intermediate	no WMS
80-0877	Advanced	no WMS
80-1583	Advanced with CELOX™	no WMS

BLEEDING CONTROL SKILLS TRAINING KITS

WITH AND WITHOUT WOUND MANAGEMENT SIMULATOR

all KITS INCLUDE:




KIT CONTENTS

- 1x Hard Storage Case with Organizer Inserts
- 20x C-A-T® Tourniquet, Trainer
- 20x 6 in. NAR® Responder ETD® (Emergency Trauma Dressing)
- 20x Just in Time Instruction Card
- 25pr Responder Nitrile Gloves, Lg
- 4x Responder Trauma Shears, 7.25 in.
- 2x NAR® ELITE™ - Red Responder Litter
- 1x **(Optional)** Wound Mgmt Simulator

+plus, by kit type:

- Basic
 - » 40x NAR® Wound Packing Gauze
- Intermediate
 - » 40x NAR® Wound Packing Gauze
 - » 20x HyFin® Chest Seal Trainer -Twin Pk
- Advanced
 - » 20x Choice of Combat Gauze® Z-Fold Hemostatic Trainer - OR - CELOX™ 5 ft. Hemostatic Trainer
 - » 20x NAR® Wound Packing Gauze
 - » 20x HyFin® Chest Seal Trainer -Twin Pk



+plus, by kit type:



BLEEDING CONTROL SKILLS TRAINING KITS

ITEM #	OPTION	WOUND MGMT SIM.
80-0880	Basic	with WMS
80-0878	Intermediate	with WMS
80-0876	Advanced	with WMS
80-1584	Advanced with CELOX™	with WMS
80-0881	Basic	no WMS
80-0879	Intermediate	no WMS
80-0877	Advanced	no WMS
80-1583	Advanced with CELOX™	no WMS

Learn more about this and other Products with a Mission® at www.NARescue.com

