

# TOMMANIKIN RIGHT ARM

## W. SHRAPNEL INJURIES

The TOMManikin® Right Arm with Shrapnel Injury features widespread penetrating shrapnel injuries, severe soft tissue damage, and bleeding from several locations. This training simulator is tourniquet responsive, this allows the users to train and utilize the Active Bleed feature of the training simulator. Constructed with 100% silicone and a durable metal bone structure, the Right Arm with Shrapnel Injuries is interchangeable with all TOMManikin® and TAMlkin® simulators.



**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



### FEATURES/SPECIFICATIONS

- Tourniquet Responsive
- Exposed tissue and muscle layers
- Integrated bleeding system
- Compatible with all TOMManikin® and TAMlkin® simulators

### **DIMENSIONS**

- L 24 in. x W 12 in. x H 6 in.
- Weight: 10 lbs.

### **TOMM RIGHT ARM - SHRAPNEL**

ITEM #	NSN
93-0085	----

\* Explore more simulator products at [www.NARescue.com](http://www.NARescue.com)