



# TOMMANIKIN RIGHT LEG

## W/ GUNSHOT WOUND

The TOMManikin® right leg with gunshot wound is tourniquet responsive and provides realistic feedback in training scenarios. Features through-and-through lateral entry wound with large medial exit wound. The 100% silicone construction with metal bone makes this training leg interchangeable with all TOMManikins® / TAMikins. Utilizing the onboard TOMManikin® bleeding system users can train on effective techniques to stop active bleeds.

**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



### FEATURES/SPECIFICATIONS

- Features through-and-through lateral entry with medial exit wound
- Tourniquet Responsive
- Compatible with all TOMManikin® and TAMikin simulators
- Wound can be used as an impaled object with a piece of wood or steel
- Available in Light Tone or Deep Tone skin

### DIMENSIONS

- Dimensions: L 36 in. x W 18 in. x D 18 in.
- Weight: 28 lbs

### TOMMANIKIN® RIGHT LEG - GSW

ITEM #	NSN
93-0097	----

\* Explore more simulator products at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277