

# TOMMANIKIN REPLACEMENT BLEEDING NECK SKINS

(10 PK)

The TOMManikin® Bleeding Neck Skins features embedded blood channels that simulate bleeding with tissue dissection. 70 operative sites within each ten pack. TOMManikin simulators provide the ability to enhance training in Tactical Combat Casualty Care (TCCC) from Point of Injury to transfer of higher care.

**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



## FEATURES/SPECIFICATIONS

- Includes 10 replacement skins
- 70 operative sites per 10 PK
- Realistic Skin-like tissue layers
- Simulated bleeding
- Works with products like the TechLine Trauma® Cric in a Box™ Simulator

## TOMMANIKIN® BLEEDING NECK SKINS (10 PACK)

ITEM #	NSN
93-0102	----

\* Explore more simulator products at [www.NARescue.com](http://www.NARescue.com)

## DIMENSIONS

- Packaged - H 10 in. x W 6 in. x L 12 in.
- Weight: 5 lbs