

TECHLINE TRAUMA RIBS IN A BOX SIMULATOR

The Techline Trauma® Wounds in a Box® series offers cost effective, compact solutions for skills training. Easy to transport, yet rugged enough for classroom or field use.

The Ribs in a Box™ task trainer is for training students in Needle Decompression of a Tension Pneumothorax. The simulated SOFTTECH®, ROCKTECH® bone and internal membrane provide realistic feedback to the trainee, while the internal cavity can accommodate an inflated exam glove or balloon to mimic escaping air when pierced. (Air Release Device not included).

Manufactured by Techline Trauma®

CONSTRUCTION

The representative skin and tissue of the Ribs in a Box™ Simulator is constructed of high-quality silicone rubber, with layers of coloration that provide a depth and realism to the product. In addition, the SOFTTECH® material has a realistic feel, and is very durable and UV resistant.

OPERATION

Contained in its own storage case, the Ribs in a Box™ is for training students in Needle Decompression of a Tension Pneumothorax. To prepare the device, the trainer places an adhesive membrane between the two layers of ROCKTECH® bones, then the assembly is pushed into the corresponding opening in the skin. If the trainer chooses to make the device express air when the needle is inserted, they may place an inflated exam glove or balloon in the cavity at its base. The device is then placed back in the storage case for use.

MAINTENANCE

The device should be rinsed in clean, warm water and allowed to dry before storage. The simulation may be removed from the storage case for cleaning. The device may be washed with mild detergents like dish soap. Allow the device to dry fully before storage to prevent mold or mildew growth.

FEATURES/SPECIFICATIONS

- Easily transportable, yet rugged enough for classroom and field use
- Patented materials give this product a realistic feel, while still being durable and UV resistant
- Simulator comes in its own storage case
- Air Release Device not included



Learn more about this and other Products with a Mission® at www.NARescue.com



DIMENSIONS

- Measurements: L 6.69 in. x W 4.63 in. x D 3.88 in.
- Weight: 1 lb. 15.8 oz.

RIBS IN A BOX™

ITEM #	NSN
93-0974	----

* Explore other Wounds in a Box® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277