

# T.A.L.O.N. MANUAL TRAINING NEEDLE SET

REV120117

## Allows for Complete Training...

Manual I/O insertion at all 7 extremity sites, including the sternum.

Includes 15 gauge, 304 stainless steel needle/catheter, sharps protector, EZ-Connect® extension set, I/O wristband, quick reference card & Sternal Locator.

Also comes with re-set tool for removal and re-use.

## Features:

- Allows for training in manual I/O insertion at all 7 extremity sites: Sternum, Bilateral proximal humerus, Bilateral proximal tibia, and Bilateral distal tibia
- Includes 15 gauge, 304 stainless steel needle/catheter, sharps protector, EZ-Connect® extension set, I/O wristband, quick reference card & Sternal Locator
- Also comes with re-set tool for removal and re-use

## Kit Contents:

- 1 x 15 gauge, 304 stainless steel needle/catheter
- 1 x Sharps Protector
- 1 x EZ-Connect® extension set
- 1 x I/O Wristband
- 1 x Quick Reference Card
- 1 x Sternal Locator
- 1 x Re-Set Tool

## Dimensions:

- Packaged: L 6 in. x W 3.25 in.
- Weight: 1.5 oz

*\*Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*



CE

## T.A.L.O.N. Manual Training Needle Set

Item #	NSN#
30-0133	

*EZ-IO Products are restricted to Federal Purchases only. Please contact Customer Service for details.*



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A  
**MISSION®**