

TOMMANIKIN® CHEST SKIN

UNINJURED

The uninjured TOMManikin® chest skin supports pleural needle decompression training scenarios. The uninjured chest skin is secured by a large section of hook and loop. This consumable item can be ordered to be bilateral chest tube compatible.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Needle decompression training skin
- No visible exterior wounds or injuries
- Compatible with all male TOMManikin® simulators
- Available in Light Tone or Deep Tone skin

UNINJURED CHEST SKIN

ITEM #	NSN
93-0148	----

** Explore more simulator products at www.NARescue.com*

DIMENSIONS

- Dimensions: L 20 in. x W 12 in. x D 6 in.
- Weight: 9 lbs