

TOMMANIKIN GUNSHOT WOUND CHEST SKIN

TOMManikin® features anterior entry GSW with non-linear (post-axillary) exit wound. Actively bleeds from both wound locations. Anterior entry wound mixes with air from onboard compressor to create bubbling bleeding effect. Supports chest seal application and pleural needle decompression. Bilateral chest tube compatible (no additional charge, requires request at time of order).



Learn more about this and other *Products with a Mission*® at www.NARescue.com



FEATURES/SPECIFICATIONS

- 100% silicone construction
- Interchangeable with all TOMManikin® simulators
- Supports Chest Seal and Pleural Needle decompression
- Non-linear (post-axillary) exit wound

TOMMANIKIN® CHEST SKIN - GSW

ITEM #	NSN
93-0149	----

* Explore more simulator products at www.NARescue.com

DIMENSIONS

- Packaged - L 20 in. x W 6 in. x H 12 in.
- Weight: 8 lbs