

# VENTIS MEDICAL HEROVENT A/C POWER SUPPLY CABLE & RECHARGEABLE BATTERY PACK

The Ventis™ Medical HeroVent® A/C Power Supply (US Cord) is a replacement power supply cord compatible with the Ventis Medical Herovent Emergency and Transport Ventilator.

## FEATURES/SPECIFICATIONS

- Replacement A/C power supply compatible with the Ventis Medical HeroVent Emergency and Transport Ventilator®
- Herovent® sold separately

## DIMENSIONS

- Dimensions: L 6 in. x W 15 in. x D 2.5 in.
- Cord Length: 47 in.
- Weight: 1.4 lbs.

The Ventis™ Medical Rechargeable Battery Pack is a replacement battery compatible with the Ventis Medical Herovent® Emergency and Transport Ventilator.

## FEATURES/SPECIFICATIONS

- **EFFICIENCY:** Settings automatically configured for common emergency situations
- **EASE OF USE:** 5 in. touchscreen interface optimized for emergency and transport care
- **DESIGN:** Turbine technology compatible with non-proprietary accessories
- **PORTABLE:** Reliable and lightweight (~3 lbs/1.5kg), with battery power for transport
- **QUALITY CARE:** Invasive & NIV modes with volume and pressure support (FiO2 and EtCO2 waveforms available)
- Compatible with the Ventis Medical HeroVent Emergency and Transport Ventilator®
- Herovent® sold separately
- Rechargeable, with a charge duration of 5+ hours
- Charging Time: 3 hours

## A/C POWER SUPPLY & BATTERY PACK

ITEM #	PRODUCT	NSN	
10-0084	A/C Power Supply	----	R <sub>X</sub> ONLY
10-0085	Battery Pack	----	R <sub>X</sub> ONLY

\* Explore other Ventis™ Medical Products at [www.NARescue.com](http://www.NARescue.com)

## DIMENSIONS

- Dimensions: L 5.2 in. x W 3.5 in. x D 1 in.
- Weight: 0.81 lb.
- Battery Information: 18V / 3.50 Ah / 63Wh



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



\*\*\* Ventis Medical HeroVent® Products are restricted to **DoW Purchases ONLY**. Please **contact Customer Service** for details.\*\*\*