

VENTIS™ MEDICAL HEROVENT ACCESSORIES

*** Ventis Medical HeroVent® Products are restricted to **DoW Purchases ONLY**. Please **contact Customer Service** for details.***

FLOW SENSOR

ITEM #
10-0090 Rx ONLY

The Ventis™ Medical HeroVent® Disposable Flow Sensor is a mass flow meter designed for proximal flow measurement in applications such as mechanical ventilation and anesthesia.

FEATURES/SPECIFICATIONS

- Replacement Flow Sensor compatible with the Ventis Medical HeroVent Emergency and Transport Ventilator®
- Herovent® sold separately
 - Disposable sensor designed for single use

DIMENSIONS

- Dimensions: L 2.73 in.
- Weight: 0.1 lbs.



FLOW SENSOR CABLE

ITEM #
10-0089 Rx ONLY

The Ventis™ Medical HeroVent® A/C Power Supply (US Cord) Flow Sensor Cable works with the HeroVent Flow Sensor for proximal flow measurement in applications such as mechanical ventilation and anesthesia.

FEATURES/SPECIFICATIONS

- Replacement Flow Sensor cable compatible with the Ventis Medical HeroVent Emergency and Transport Ventilator®
- Herovent® sold separately

DIMENSIONS

- Cord Length: 83 in.
- Weight: 0.17 lbs.



ETCO2 SENSOR

ITEM #
10-0088 Rx ONLY

The Ventis™ Medical HeroVent® EtCO2 sensor provides real-time monitoring of end-tidal carbon dioxide, or EtCO2, in order to evaluate ventilation and track patient response.

FEATURES/SPECIFICATIONS

- Replacement EtCO2 sensor compatible with the Ventis Medical HeroVent Emergency and Transport Ventilator®
- Herovent® sold separately

DIMENSIONS

- Weight: 0.11 lbs.



Learn more about this and other Products with a Mission® at www.NARescue.com



VENTIS™ MEDICAL HEROVENT ACCESSORIES

DoW PURCHASES ONLY

OXYGEN ENRICH. KIT

ITEM #	
10-0093	R _X ONLY

The Ventis™ Medical HeroVent® Oxygen Enrichment Kit is a reusable replacement kit compatible with the Ventis Medical Herovent Emergency and Transport Ventilator.

FEATURES/SPECIFICATIONS

- Reusable kit designed for multiple uses
- Herovent® sold separately

DIMENSIONS

- Dimensions: L 2.73 in.
- Weight: 0.1 lbs.



BREATH EXHALATION VALVE

ITEM #	
10-0094	R _X ONLY

The Ventis™ Medical HeroVent® Breath Exhalation Valve 22mm with Monitoring Line - 1.8m is a replacement valve compatible with the Ventis Medical Herovent Emergency and Transport Ventilator.

FEATURES/SPECIFICATIONS

- Replacement exhalation valve compatible with the Ventis Medical HeroVent® Emergency and Transport Ventilator®
- Herovent® sold separately

DIMENSIONS

- Rolled Circumference: 6 in.
- Length of Cable: 1.8m (5.91 ft.)
- Weight: 0.2 lbs.



INLET FILTER 5-PK

ITEM #	
10-0091	R _X ONLY

The Ventis™ Medical HeroVent® Inlet Filters are replacement filters compatible with the Ventis Medical Herovent Emergency and Transport Ventilator.

FEATURES/SPECIFICATIONS

- Herovent® sold separately
- Filters come in a pack of 5

DIMENSIONS

- Circumference (per 1 Filter): 115 in.
- Weight (per 1 Filter): 0.09 g.
- Weight (per 5pk): 0.48 g.



Learn more about this and other Products with a Mission® at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277