

SMRT™ TASK TRAINER - SURGICAL AIRWAY FRONT

STRATEGIC MEDICAL RESEARCH & TRAINING™

The SMRT Surgical Airway Front is a replacement skin designed to work with the SMRT Surgical Airway Moulage Task Trainer.

SMRT surgical airway (Cric) trainer is designed to make cricothyroidotomy training realistic, fast and easy. Whether you're are training Physicians, Paramedics or Military Medics, this device will meet your "cric" training needs. Resetting this replacement skin takes seconds making it perfect for skills stations or scenario training. Complete with realistic textured sacrificial skin and blood packs that give a realistic feel and response.

This flexible collar has integrated stainless steel mesh to protect the wearer during the open surgical airway procedure. This device is not designed for needle or wire-guided procedures. This device should be used in an observed setting with trained health professionals familiar with performing an open cricothyroidotomy.

WARNING: This device is designed to be used only for the open surgical airway procedure. No needle or wire-guided procedures should be performed on this training device. A trained medical professional should observe the trainee in a 1:1 ratio to ensure the safety and accuracy of the procedure.

The innovative and wearable SMRT™ Moulage Task Trainers from Strategic Medical Research & Training™ create an unmatched experience of realism featuring realistic moulage and active bleeding across a wide spectrum of woundings and interventions. Made with durable materials and designed for a "one size fits most" range of limb and body sizes, the SMRT™ Moulage Task Trainers are engineered to work with live role players as well as to augment existing simulation systems—while allowing flexibility in placement and utilization. All SMRT™ Moulage Task Trainers work with SMRT™ Trauma Tables, SMRT™ Bleeding Systems, and common hydration systems, making them standalone and ready to use, while also allowing ease of integration into existing training systems. Utilized by Military, Law Enforcement, EMS, and Hospitals to create unsurpassed realism in hemorrhage control training, SMRT™ Training Aids are realistic, cost effective, highly mobile, durable, and low maintenance.

FEATURES/SPECIFICATIONS

- Surgical Airway Front is the replacement skin for the SMRT™ Surgical Airway Moulage Task Trainer
- Available in Light Tone or Deep Tone
- Unsurpassed realism in hemorrhage control training
- One size fits most ranges of body sizes

Learn more about this and other Products with a Mission® at www.NARescue.com



Available in light and deep tone.

SMRT™ SURGICAL AIRWAY FRONT

ITEM #	COLOR
93-1153	Light Tone
93-1154	Deep Tone

DIMENSIONS

- L 13 in. x W 6 in. x D 0.25 in.
- Weight: 0.3 lbs.

* Explore other SMRT™ moulage task trainers at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277