

SMRT™ TASK TRAINER - TRAUMA TABLE - HYPERLITE

STRATEGIC MEDICAL RESEARCH & TRAINING™

The SMRT™ Hyperlite Trauma Table is the solution to your hemorrhage control training needs. This system allows for unsurpassed realistic, high-pressure simulated bleeding that is compatible with the full suite of SMRT-simulated wounds.


This patented system collects, filters, warms, and delivers high-pressure simulated blood to a full-sized manikin. Utilize your existing manikin by applying SMRT's role-player wounds to create an ultra-realistic, high fidelity hemorrhage control task trainer.

The Hyperlite configuration allows for ease of storage and shipping making this the ideal system for mobile training sites as well as simulation labs.

This model utilizes a folding table and flexible cover system.

The innovative and wearable SMRT™ Moulage Task Trainers from Strategic Medical Research & Training™ create an unmatched experience of realism featuring realistic moulage and active bleeding across a wide spectrum of woundings and interventions. Made with durable materials and designed for a "one size fits most" range of limb and body sizes, the SMRT™ Moulage Task Trainers are engineered to work with live role players as well as to augment existing simulation systems—while allowing flexibility in placement and utilization. All SMRT™ Moulage Task Trainers work with SMRT™ Trauma Tables, SMRT™ Bleeding Systems, and common hydration systems, making them standalone and ready to use, while also allowing ease of integration into existing training systems. Utilized by Military, Law Enforcement, EMS, and Hospitals to create unsurpassed realism in hemorrhage control training, SMRT™ Training Aids are realistic, cost effective, highly mobile, durable, and low maintenance.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Unsurpassed realism in hemorrhage control training
- Trauma Table - Hyperlite Variant is compatible with other SMRT™ Task Trainers
- Up to eight (8) bleeds at one time
- Recirculates and warms simulated blood for realism

SMRT™ TRAUMA TABLE - HYPERLITE

ITEM #

93-1212

* Explore other SMRT™ moulage task trainers at www.NARescue.com

DIMENSIONS

- Packaged Dimensions: L 36 in. x W 36 in. x D 20 in.
- Weight: 121.6 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

SMRT™ TASK TRAINER - TRAUMA TABLE - HYPERLITE

STRATEGIC MEDICAL RESEARCH & TRAINING™

