

MODULAR TRAINING SYSTEM

Module One: Skill Stations

REV040218

Train up to 20 Students at Individual Skill Stations...

The TCCC/TECC Modular Training System from North American Rescue provides all the necessary tactical medical and casualty extraction/evacuation equipment needed when teaching the PHTLS TCCC Certification required tasks for completing the Tactical Combat Casualty Care (TCCC) course. This system is available as a turnkey solution for training centers and educational institutions. Each training kit is designed to equip up to twenty (20) students and consists of three modules:

Module One – Skill Stations: Provides the necessary equipment to perform PHTLS TCCC Certification required tasks when training at individual skill stations. Includes all the necessary equipment to instruct up to twenty (20) students in the treatment of the leading causes of potentially preventable death in the tactical environment. Contents are packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory. Addresses extremity hemorrhage, wound packing, chest seal application, needle chest decompression, junctional hemorrhage, basic airway adjuncts, surgical airway, splinting, hypothermia prevention & management, IV skills and eye injuries.

Features:

- Provides the necessary equipment to train up to twenty (20) students at individual skill stations in the treatment of the leading causes of potentially preventable death in the tactical environment
- Packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory (Resupply kit with contents only is also available)
- Addresses extremity hemorrhage, wound packing, chest seal application, needle chest decompression, junctional hemorrhage, basic airway adjuncts, surgical airway, splinting, hypothermia prevention & management, IV skills and eye injuries

Dimensions:

- Hard Storage Case:
External – L 30.75 in. x W 20.5 in. x D 11.6 in.
Internal – L 28.5 in. x W 17.5 in. x D 10.5 in.
- Weights: Module One - 53 lb

TCCC Modular Training System - Module One

Module Options	Item #	NSN#
Contents in Rigid Case	85-0519	
Re-Supply (Contents Only)	85-1900	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
— WITH A —
MISSION®

MODULAR TRAINING SYSTEM *Module One: Skill Stations*

REV040218



Module One Contents:

- 20 x C-A-T® Tourniquet, Trainer
- 20 x NAR 6 in. Emergency Trauma Dressing
- 20 x NAR S-Rolled Gauze
- 20 x Combat Gauze Trainer
- 1 x JETT® Junctional Emergency Treatment Tool
- 5 x Nasopharyngeal Airway, 28F
- 1 x Cyclone® BVM
- 2 x Tactical CricKit®
- 10 x HyFin® Twin Pack Dressing, Trainer
- 5 x ARS® Needle Decompression Kit, 14 ga.
- 20 x Saline Lock Kit
- 5 x NAR Sharps Container
- 1 x HPMK® Hypothermia Prevention and Management Kit
- 5 x SAM® Splint
- 5 x NAR Triangular bandage
- 5 x Polycarbonate Eye Shield
- 1 x Hard Storage Case

Hard Case Specifications:

- Retractable extension handle
- 4 strong polyurethane wheels with stainless steel bearings
- Easy-open double throw latches
- Open cell core with solid wall design – strong & lightweight
- O-ring seal
- Automatic pressure equalization valve
- Fold down handles
- Stainless steel hardware and padlock protectors
- Desert tan color



ReSupply Kit with contents only is also available



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®