

CELOX RAPID HEMOSTATIC GAUZE

Tactical Medicine means working in urgent and possibly life-threatening conditions. Celox™ RAPID bridges the significant capability gap between this urgency and the current standard of treatment for hemorrhage control. Celox RAPID offers a step change in speed, a breakthrough that expands the utility of hemostatics in Tactical Field Care.

Gauze contains an x-ray detectable line embedded within the material allowing for identification of the gauze under standard imaging to ensure quick and easy removal from the patient.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Includes embedded X-Ray detectable line
- Stops severe arterial bleeding with only one minute of compression applied
- Compact Z-folded format making it easy to use & carry
- Generates no heat

CELOX RAPID GAUZE

ITEM #	NSN	
30-0069	6510-01-623-9903	Rx ONLY

* Explore more gauze options at www.NARescue.com

DIMENSIONS

- Packaged: H 5.25 in. x W 4.25 in. x D 0.5 in.
- Gauze: W 3 in. x 5 ft
- Weight: 11 oz