

# K9 MINI MEDIC KIT

**Your K-9 partner would defend you with his life; don't leave his to chance.**

The K9 Mini Medic Kit provides a compact, versatile platform equipped with the essential medical equipment to treat multiple traumatic injuries, such as penetrating and blast-type wounds, which are common to K9s and their handlers in their operational arenas. The K9 specific gear enables handlers to render care to an injured canine when veterinary care is not immediately available.

This optimally sized Mini Medic bag features external pouches and side/ bottom MOLLE to allow for external tourniquet or other attachment options. A large pouch on back panel includes a drop-down pouch with abrasive resistant padded bottom providing extra storage for more gear such as compact breaching tools. It's external expansion zipper allows for increased load capacity if needed. The inside of the bag has two small and one large, clear tear-away pouches for medical or other equipment supply storage. Along with two mesh inner pouches and an internal pouch on the inner lid, many storage options exist.

## FEATURES/SPECIFICATIONS

- Compact, multi-purpose trauma aid bag with life-saving trauma equipment and K9 specific gear for point of wounding treatment of multiple traumatic injuries
- K9 specific gear enables handlers to render care to an injured canine when veterinary care is not immediately available
- Versatile platform features two external pockets with two external tourniquet holders, one external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack and an expansion zipper to increase load capacity
- Additional configuration capabilities include external MOLLE on sides and bottom of bag, internal MOLLE rows on inside back of bag, three clear tear away pouches / two internal mesh pouches on the inside main compartment, and one larger internal pouch on inner lid
- Internal tape roll strap also included



## DIMENSIONS

- H 14 in. x W 11 in. x D 5 in. (Depth expandable to 7.25 in.)
- Weight: 12.5 lb

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## MINI MEDIC KITS

COLOR	LEVEL	ITEM #	NSN #
Black	Basic	80-1448	---
Coyote	Basic	80-1449	---
OD Green	Basic	80-1451	---
RGR Green	Basic	80-1452	---
Multicam	Basic	80-1450	---
Black	Combat Gauze	80-1453	---
Coyote	Combat Gauze	80-1454	---
OD Green	Combat Gauze	80-1456	---
RGR Green	Combat Gauze	80-1457	---
Multicam	Combat Gauze	80-1455	---
Black	CELOX™	80-1463	---
Coyote	CELOX™	80-1464	---
OD Green	CELOX™	80-1466	---
RGR Green	CELOX™	80-1467	---
Multicam	CELOX™	80-1465	---
Black	Adv. w/ Combat Gauze	80-1458	---
Coyote	Adv. w/ Combat Gauze	80-1459	---
OD Green	Adv. w/ Combat Gauze	80-1461	---
RGR Green	Adv. w/ Combat Gauze	80-1462	---
Multicam	Adv. w/ Combat Gauze	80-1460	---
Black	Adv. w/ CELOX™	80-1468	---
Coyote	Adv. w/ CELOX™	80-1469	---
OD Green	Adv. w/ CELOX™	80-1471	---
RGR Green	Adv. w/ CELOX™	80-1472	---
Multicam	Adv. w/ CELOX™	80-1470	---

Rx ONLY

# K9 MINI MEDIC KIT

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



Every K9 MINI MEDIC KIT INCLUDES THESE CORE ITEMS:

 QTY 2 C-A-T® TQ Holder	 QTY 4 C-A-T® Tourniquet	 QTY 2 NAR® Gauze S-Rolled	 QTY 4 NAR® Mini Responder 4" ETD®	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 4 Petrolatum Gauze
 QTY 8 NAR Alcohol Prep Pads	 QTY 10 Flexible Fabric Bandages	 QTY 1 2" NAR® Surgical Tape	 QTY 1 Emergency Survival Wrap	 QTY 1 Tactical Muzzle	 QTY 1 NAR® K-9 Litter
 QTY 4 Green ChemLight™	 QTY 1 Infrared Square	 QTY 1 Streamlight Headlamp	 QTY 1 NAR® Trauma Shears, Lg	 QTY 1 Veterinary Thermometer	
 PAIR 5 Responder Black Gloves, LG	 QTY 1 Permanent Marker, Large	 QTY 1 Patch 2" x 2" Medical Cross			

+plus, by kit type:

<b>Basic</b>  QTY 2 NAR® Gauze S-Rolled	- Or -	<b>Combat Gauze</b>  QTY 2 Combat Gauze® Z-Fold	- Or -	<b>CELOX™</b>  QTY 2 CELOX™ Gauze Z-Fold
---	--------	---	--------	--

<b>Advanced</b>  QTY 4 Choice of Hemostatic	 QTY 2 Pre-lubed NAR® Naso, 28F	 QTY 2 ARS® Needle, 14g
---	--	----------------------------------



Designed to provide various capabilities for different mission profiles...