

BUTTERFLYBVM™

The ButterflyBVM™ limits tidal volumes and helps medical personnel give air according to published guidelines.

The ButterflyBVM comes in three options: **Standard** (Adult mask only), **Plus** (Pediatric & Adult masks with PEEP/Filter and elbow), or **Ready** (Neonatal, Infant, Pediatric, & Adult masks with PEEP/Filter and elbow)

Learn more about this and other *Products with a Mission*® at www.NARescue.com



FEATURES/SPECIFICATIONS

Adjustable Tidal Volume

- Giving patients a precise volume of air is not easy, but with a turn of a wheel you can restrict the ButterflyBVM's maximum air output to match your patient's needs. Follows Broselow and Handtevy groups and comes standard on every ButterflyBVM model
- By restricting tidal volumes to safer levels, the risks to patients of hyperventilation are reduced

Variable Pop-Off

- Worried about barotrauma? Our PIP dial gives lifesavers unprecedented control over the pressures that they deliver to patients without the need for a manometer

Rate Control

- Have you ever seen someone squeeze a BVM too rapidly? The ButterflyBVM's built-in rate control feature discourages hyperventilation by excessive rate.

Accessory Ready

- The ButterflyBVM is ready to adapt to complex scenarios by receiving standard PEEP valves, filters, and ET-CO2 detectors while also fitting all standard patient interfaces

Internal O2 Reservoir

- With the help of its internal O2 reservoir, the ButterflyBVM delivers a 90%+ FIO2 to patients of all sizes when supplemental O2 is introduced at 15LPM



BUTTERFLYBVM™ Rx ONLY

ITEM #	OPTION
10-0099	Standard
10-0100	Plus
10-0101	Ready

* Explore other BVMs and respiratory products at www.NARescue.com

DIMENSIONS

- Packaged Dimensions: L 17 in. x W 12 in. x D 6 in.

STANDARD:

- Packaged Weight: 1.2 lb.

PLUS:

- Packaged Weight: 1.33 lb.

READY:

- Packaged Weight: 1.49 lb.

BUTTERFLYBVM™

BUTTERFLYBVM™ STANDARD:

- Tidal Volume Control
- Peak Inspiratory Pressure Control
- Rate Control
- Adult-Peds-Infant BVM in one
- Oxygen Tubing
- Adult Mask

BUTTERFLYBVM PLUS™:

- Tidal Volume Control
- Peak Inspiratory Pressure Control
- Rate Control
- Adult-Peds-Infant BVM in one
- Oxygen Tubing
- Adult Mask
- Pediatric Mask
- PEEP Valve
- Filter
- Ergonomic Elbow

BUTTERFLYBVM READY™:

- Tidal Volume Control
- Peak Inspiratory Pressure Control
- Rate Control
- Adult-Peds-Infant BVM in one
- Oxygen Tubing
- Adult Mask
- Pediatric Mask
- PEEP Valve
- Filter
- Ergonomic Elbow
- Infant Mask
- Neonatal Mask



Learn more about this and other Products with a Mission® at www.NARescue.com



BUTTERFLYBVM™

ITEM #	OPTION
10-0099	Standard
10-0100	Plus
10-0101	Ready

* Explore other BVMs and respiratory products at www.NARescue.com