

MODULAR TRAINING SYSTEM

Module Two: IFAKs

20 Pre-packed TCCC Training IFAKs...

The TCCC/TECC Modular Training System from North American Rescue provides all the necessary tactical medical and casualty extraction/evacuation equipment needed when teaching the PHTLS TCCC Certification required tasks for completing the Tactical Combat Casualty Care (TCCC) course. This system is available as a turnkey solution for training centers and educational institutions. Each training kit is designed to equip up to twenty (20) students and consists of three modules:

Module Two – IFAKs: Provides twenty training IFAKs for student use during trauma lanes or scenarios. Each IFAK contains the necessary equipment for initial TCCC/TECC care in an easy to attach modular carrying bag that can be worn as a drop-down leg rig attached to a belt with a quick-detach buckle system. The IFAK also has a removable belt and leg strap so it can be MOLLE attached to a vest as another carrying option. The kit contains 20 pre-packed TCCC Training IFAKs in a hard, wheeled case that allows for easy transport and storage.

Features:

- Contains 20 pre-packed TCCC Training IFAKs for use by students during trauma lanes or training scenarios.
- Individual First Aid Kits (also available for separate purchase) come packaged in a self contained wheeled hard, case that allows for easy transport and storage.

Dimensions:

- Hard Storage Case:
External – L 30.75 in. x W 20.5 in. x D 11.6 in.
Internal – L 28.5 in. x W 17.5 in. x D 10.5 in.
- Weights: Module Two - 59 lb



TCCC Modular Training System - Module Two

Item #

NSN#

85-0520

*Twenty training IFAKs for student use
during trauma lanes or scenarios.*

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

Module Two Contents:

- 1x Hard Storage Case
- 20x TCCC IFAKS

Each IFAK contains the following:


- 1x C-A-T® Tourniquet, Trainer
- 1x NAR 6 in. Emergency Trauma Dressing
- 1x NAR S-Rolled Gauze
- 1x Combat Gauze Trainer
- 1x HyFin® Twin Pack Dressing, Trainer
- 1x Nasopharyngeal Airway, 28f
- 1x ARS®, Needle Decom Kit, 14 ga.

- 1x Permanent Marker
- 1x Combat Casualty Card
- 2x Bear Claw™ Glove Kit, Large
- 1x 2 in. Surgical Tape



Hard Case Specifications:

- Retractable extension handle
- 4 strong polyurethane wheels with stainless steel bearings
- Easy-open double throw latches
- Open cell core with solid wall design – strong & lightweight
- O-ring seal
- Automatic pressure equalization valve
- Fold down handles
- Stainless steel hardware and padlock protectors
- Desert tan color



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277