

MODULAR TRAINING SYSTEM

Module Three: Extraction

A Turn-Key Solution for Training Centers...

The TCCC/TECC Modular Training System from North American Rescue provides all the necessary tactical medical and casualty extraction/evacuation equipment needed when teaching the PHTLS TCCC Certification required tasks for completing the Tactical Combat Casualty Care (TCCC) course. This system is available as a turnkey solution for training centers and educational institutions. Each training kit is designed to equip up to twenty (20) students and consists of three modules:

Module Three – Extraction: Provides various high-threat extraction devices and casualty evacuation systems to train students in different options to move casualties during the all three TCCC Phases of Care. Kit includes High Threat Extraction drag devices, versatile casualty movement Hasty Harness, poleless/non-rigid litter option and the collapsible quad-fold rigid evacuation platform standard – Talon II 90C Litter.

Features:

- Provides various high threat extraction devices and casualty evacuation systems to train students in different options to move casualties during the all three TCCC Phases of Care.
- Packaged in a self-contained, wheeled, rigid case that allows for easy transport and storage.

Dimensions:

- Hard Storage Case:
External – L 30.75 in. x W 20.5 in. x D 11.6 in.
Internal – L 28.5 in. x W 17.5 in. x D 10.5 in.
- Weights: Module Three – 45 lb



TCCC Modular Training System - Module Three

Item #

NSN#

85-0521

Provides all the necessary tactical medical and casualty extraction/evacuation equipment needed when teaching the PHTLS TCCC Certification

PRODUCTS WITH A
MISSION[®]



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

Module Three Contents:

- 1x Hard Storage Case
- 1x Talon II® 90C Collapsible Handle Litter
- 1x EEL, Emergency Evacuation Litter- BLK
- 1x TED, Tactical Extraction Device- BLK
- 1x Hasty Harness - Blk
- 1x RAT Strap, Rescue Assault Tether
- 1x Dragon Handle System
- 1x Shield Carabiner



Hard Case Specifications:

- Retractable extension handle
- 4 strong polyurethane wheels with stainless steel bearings
- Easy-open double throw latches
- Open cell core with solid wall design – strong & lightweight
- O-ring seal
- Automatic pressure equalization valve
- Fold down handles
- Stainless steel hardware and padlock protectors
- Desert tan color



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277