

SMRT™ TASK TRAINER - ANIMAL ATTACK

STRATEGIC MEDICAL RESEARCH & TRAINING™

The SMRT™ Animal Attack Moulage Task Trainer is designed to simulate a wild animal attack wound or injury from a police k9 or military working dog. The task trainer can be placed on the upper arm, thigh, and forearm in adults or thighs on pediatric role players.

The innovative and wearable SMRT™ Moulage Task Trainers from Strategic Medical Research & Training™ create an unmatched experience of realism featuring realistic moulage and active bleeding across a wide spectrum of woundings and interventions. Made with durable materials and designed for a "one size fits most" range of limb and body sizes, the SMRT™ Moulage Task Trainers are engineered to work with live role players as well as to augment existing simulation systems—while allowing flexibility in placement and utilization. All SMRT™ Moulage Task Trainers work with SMRT™ Trauma Tables, SMRT™ Bleeding Systems, and common hydration systems, making them standalone and ready to use, while also allowing ease of integration into existing training systems. Utilized by Military, Law Enforcement, EMS, and Hospitals to create unsurpassed realism in hemorrhage control training, SMRT™ Training Aids are realistic, cost effective, highly mobile, durable, and low maintenance.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Unsurpassed realism in hemorrhage control training
- One size fits most range of limb and body sizes
- Animal Attack compatible with other SMRT™ Task Trainers
- Available in Light Tone or Deep Tone

SMRT™ ANIMAL ATTACK

ITEM #	COLOR
93-1201	 Light Tone
93-1202	 Deep Tone

* Explore other SMRT™ moulage task trainers at www.NAResce.com

DIMENSIONS

- L 9 in. x W 5 in. x D 0.5 in.
- Weight: 0.5 lbs.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277