

TALON® BACKREST

Compatible with the Talon®, Raven®, and MedEVAC litters

For appropriate treatment and improved comfort, supine is not always the best position for patients. When casualties can be better served by being placed at a 45° or 90° angle, the Talon® Series Backrest from North American Rescue® offers an efficient and reliable solution.

This litter backrest attaches to the 1.5 in. poles of the Talon® system, the premiere litter system on the market. The backrest is also engineered to fit the Raven® Model 90C as well as the Medevac Litter.

The backrest cover is made of a Nomex/Kevlar blend: Nomex for flame-resistance and heat-energy absorption; Kevlar for abrasion-resistance and long life.

Learn more about this and other **Products with a Mission®** at www.NARescue.com



NOTE: All accessories are sold separately. Litters are not included.



When the patient's condition requires a 45° or 90° angle, rely on the Talon® Series Backrest.

DIMENSIONS

- H 31.5 in. x W 22.5 in. x D 1.6 in.
- Weight: 2 lb 13 oz

FEATURES/SPECIFICATIONS

- Allows patients to be placed at 90° or 45° angles
- Made of Nomex/Kevlar-blend fabric
- Compatible for use with the Talon®, Raven® and MedEVAC litters
- Firm, comfortable, and secure

TALON® BACKREST

ITEM #	NSN
60-0027	6530-01-594-7444
60-0027-EU	6530-01-594-7444



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277