

OSL™ BACKREST

Learn more about this and other **Products with a Mission®** at www.NARescue.com



The OSL Backrest attaches to the 2 in. poles of the OSL litter system allowing the patient to be placed at both a 90° and a 45° angle. The OSL backrest cover material is a very durable Nomex/ Kevlar blend.



NOTE: All accessories are sold separately. Litters are not included.

FEATURES/SPECIFICATIONS

- SEL#: 09ME-01-COTS
- Covered with durable Nomex/ Kevlar-blend fabric
- Positions patients at 45° and 90° angles
- Attaches to the 2 in. poles of the OSL™ Litter
- Increases patient comfort

OSL™ BACKREST

ITEM #	NSN
60-0024	NATO NSN AVAIL
60-0024-EU	----



DIMENSIONS

- H 31.5 in. x W 30.25 in. x D 1.6 in.
- Weight: 2 lb 15 oz

