

# MINI TRAUMA SHEET KIT

## Easy Access...

Introducing the NAR® Mini Trauma Sheet Kit, a portable, compact rollable sheet packed with essential first aid and essential life-saving equipment that addresses injuries caused by penetrating or blast trauma. The Mini Trauma Sheet Kit can be rolled up to be carried to an incident, or hung in an ambulance, vehicle, or aircraft for easy access to medical supplies.

The Mini Trauma Sheet Kit features clear vinyl pouches for easy identification of contents, and each pouch can be accessed through the large zipper opening. Pouches have exterior loop fields enabling the use of hook and loop name tape labeling. The kit includes two "BLEED", one "MISC/MEDS", and one "AIRWAY" name tapes.

The Mini Trauma Sheet Kits are available in different levels of care: Basic, and Basic with QuikClot® Bleeding Control Dressing. The Kits are available in your choice of the following nylon colors: Black, Red, and Coyote.



## FEATURES/SPECIFICATIONS

- Can be rolled up and carried with the attached handles, or opened and hung using the sewn in D-Rings
- Four large pouches and one lower profile zippered pocket for maximum capacity

- Compact design for smaller vehicles, ambulances, etc.
- Ability to pack out however you like for best possible outcomes
- Hook and Loop attachment for name tape additions for better quality quick visual access



Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



## MINI TRAUMA SHEET KIT

ITEM #	COLOR	OPTION	NSN
80-1808	Black	Basic	----
80-1809	Coyote	Basic	----
80-1807	Red	Basic	----
80-1805	Black	Basic with BCD	----
80-1806	Coyote	Basic with BCD	----
80-1804	Red	Basic with BCD	----


## DIMENSIONS

- Closed: W 13 in. x L 9.75 in. x D 5 in.
- Open: W 13 in. x L 19 in. x D 2 in.
- Weight: 2.13 lbs.

# MINI TRAUMA SHEET KIT



**+plus:**



\* RED Mini Trauma Sheet kits are packed with Orange C-A-T® Tourniquets and RED/BLK name patches.

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## KIT CONTENTS

**+plus**  
w/ BCD

- 2 x C-A-T® Tourniquet - Orange
- 2 x NAR® Responder ETD® - 4 in. (Emergency Trauma Dressing)
- 2 x QuikClot® Bleeding Control Dressing - 3 in. x 4 yd. Z-Fold
- 2 x NAR® Compressed Gauze - 4.5 in. x 4.1 yd.
- 2 x NAR® Nasopharyngeal Airway - 28FR - Pre-lubricated
- 2 x HyFin® Vent Chest Seal Dressing - Twin Pack
- 2 x Gauze Pad (2 pack) - 4 in. x 4 in.
- 1 x NAR® Flat Tape - 2 in.
- 1 x Elastic Wrap Bandage - 4 in.
- 1 x Minor Care Bandages Kit
  - » 15 x Bandage - 1 in. x 3 in.
  - » 5 x Knuckle Bandage - 1.5 in. x 3 in.
  - » 5 x Flexible Finger Tip Bandage
  - » 5 x Flexible Fabric Bandages
  - » 1 x Wound Closure Strip (10 pack)
- 1 x NAR® Survival Blanket
- 2 x CPR Microshield Mask
- 4 x Sting & Bite Antiseptic
- 4 x Antibiotic Ointment
- 1 x NAR® Mini Duct Tape - 2 in.
- 1 x NAR® Trauma Shears - Small
- 1 x Permanent Marker - Large
- 2pr Nitrile Gloves - Large
- 1 x Name Patch Set
  - » 2x Bleed Name Patch
  - » 1x Misc / Meds Name Patch
  - » 1x Airway Name Patch