

# MINI TRAUMA SHEET KIT - CLASS A

**ANSI Standards for workplace first aid kits have been updated! DO NOT FALL OUT OF COMPLIANCE!**

**Updated Items Included In Class A Standards:**

- 1 x HyFin® Vent Chest Seal Twin Pack
- 2 x Sterile Oval Eye Pad
- 10 x Hand Sanitizer, 3g
- 1 x NAR® Survival Blanket - 56 in. x 84 in.
- 1 x Black Permanent Marker, Large

The Mini Trauma Sheet Kit - Class A features clear vinyl pouches for easy identification of contents, and each pouch can be accessed through the large zipper opening. The Mini Trauma Sheet Kit - Class A from North American Rescue® is designed as a portable First Aid Kit for the workplace that exceeds ANSI/ISEA Z308.1-2021 Standard for Class A kits. Our kits exceed the ANSI and OSHA standards by adding additional bleeding control elements. Packaged in a convenient rolled up sheet that is easily transported to the patient's side, this kit offers the ability to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T® Tourniquet, the fastest, safest, most effective prehospital field tourniquet in the world.



**FEATURES/SPECIFICATIONS**

- Workplace first aid kits that exceed ANSI/ISEA Z308.1-2021 Standard for Class A kits
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T® Tourniquet, the fastest, safest, most effective prehospital field tourniquet in the world
- Designed so the kit can be easily carried to the patient's side
- Can be rolled up and carried with the attached handles, or opened and hung using the sewn in D-Rings
- Four large pouches and one lower profile zippered pocket for maximum capacity
- Compact design for vehicles

**DIMENSIONS**

- Closed: W 13 in. x L 9.75 in. x D 5 in.
- Open: W 13 in. x L 19 in. x D 2 in.
- Weight: 2 lbs. 2 oz.

**MINI TRAUMA SHEET KIT - CLASS A**

| ITEM #  | COLOR  | OPTION      | NSN  |
|---------|--------|-------------|------|
| 80-1926 | Black  | Basic       | ---- |
| 80-1929 | Coyote | Basic       | ---- |
| 80-1911 | Red    | Basic       | ---- |
| 80-1927 | Black  | with BCD    | ---- |
| 80-1930 | Coyote | with BCD    | ---- |
| 80-1912 | Red    | with BCD    | ---- |
| 80-1928 | Black  | with CELOX™ | ---- |
| 80-1931 | Coyote | with CELOX™ | ---- |
| 80-1913 | Red    | with CELOX™ | ---- |


**Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)**



# MINI TRAUMA SHEET KIT - CLASS A



plus by kit type:



## KIT CONTENTS

- 1x C-A-T® Tourniquet - Orange
- 1x NAR® Compressed Gauze - 4.5 in. x 4.1 yd.
- 1x NAR® Flat Responder Emergency Trauma Dressing (ETD®) - 4 in.
- 1x HyFin® Vented Chest Seal Dressing (Twin Pack)
- 1x NAR® Survival Blanket - 56 in. x 84 in.
- 1x Triangular Bandage - 40 in. x 40 in. x 56 in.
- 1x SAM® Splint
- 2x Gauze Pad (2 pack) - 4 in. x 4 in.
- 2x Gauze Pad - 5 in. x 9 in.
- 1x BurnTec® Dressing - 4 in. x 4 in.
- 1x CPR Microshield Mask
- 10x Burn Gel
- 10x Antiseptic Towelette
- 10x Antiseptic Ointment
- 10x Hand Sanitizer
- 1x Instant Cold Pack - 5 in. x 5.5 in.
- 1x Eye Wash Solution - 1 oz.
- 2x Sterile Oval Eye Pad
- 1x Elastic Wrap Bandage - 2 in. x 5 yd.
- 1x Elastic Wrap Bandage - 4 in. x 5 yd.
- 1x Surgical Tape - 1 in. x 10 yd.
- 16x Fabric Bandages - 1 in. x 3 in.
- 2 pr Responder Gloves - Large
- 1 pr Responder Trauma Shears - 7.25 in.
- 1x Permanent Marker - Large
- 1x First Aid Pocket Guide
- 1x Just In Time Bleeding Control Instructions
- 1x Trauma & First Aid Kit Contents Inventory Card
- 1x Inspection Card

plus by kit type:

- 1x NAR® Compressed Gauze - 4.5 in. x 4.1 yd. - OR -
- 1x QuikClot® Bleeding Control Dressing - 3 in. x 4 yd. - OR -
- 1x Celox™ First Aid Z-Fold Hemostatic Dressing - 3 in. x 5 yd.

## MINI TRAUMA SHEET KIT - CLASS A

| ITEM #  | COLOR  | OPTION      | NSN  |
|---------|--------|-------------|------|
| 80-1926 | Black  | Basic       | ---- |
| 80-1929 | Coyote | Basic       | ---- |
| 80-1911 | Red    | Basic       | ---- |
| 80-1927 | Black  | with BCD    | ---- |
| 80-1930 | Coyote | with BCD    | ---- |
| 80-1912 | Red    | with BCD    | ---- |
| 80-1928 | Black  | with CELOX™ | ---- |
| 80-1931 | Coyote | with CELOX™ | ---- |
| 80-1913 | Red    | with CELOX™ | ---- |

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

