

EZ-TOURNIQUET POUCH (GEN 2)

Flawless Deployment for Critical Moments

Stay mission-ready with the EZ-Tourniquet Pouch (Gen. 2), your ultimate one-pull solution for life-saving response. Built from premium Cordura 500D/500D and the highest-grade materials, this new generation pouch is designed to withstand the harshest conditions, ensuring reliability when it matters most.

Redesigned for enhanced versatility, the Gen. 2 seamlessly integrates with a wide range of tactical gear through dual Velcro straps or the MOLLE-compatible strap. This adaptable design keeps your tourniquet always within reach, whether attached to belts, vests, or other setups. Trusted by elite forces worldwide, its rugged durability and rapid deployment capability make it an indispensable tool in high-stakes operations. Be prepared—every second counts.

Seamless Integration, Life-Saving Speed

The development of the EZ-Tourniquet Pouch was a concerted effort, informed by the firsthand experiences of military experts and field operatives actively engaged in high-stakes environments. Working closely with seasoned combat veterans and emergency medicine specialists, we refined every aspect of the pouch's design to meet the precise demands of real-world operations. This collaboration has ensured that the pouch not only meets but exceeds the rigorous expectations of those who face critical situations daily, providing a tool that is both reliable and indispensable in saving lives.

FEATURES/SPECIFICATIONS

- One-Pull Activation:**
The patented design enables immediate tourniquet deployment with a single pull. Even in situations where one hand is incapacitated, the tourniquet can be effortlessly deployed within seconds—ensuring swift, life-saving action in critical moments.
- Combat Durability:**
Built from top-tier Cordura 500D/500D and advanced military-grade materials, the Gen. 2 pouch offers exceptional resilience. It withstands the harshest environments, from arid deserts to extreme humidity, guaranteeing peak performance through countless missions.
- Enhanced Versatility:**
Now featuring Velcro straps, the EZ-Tourniquet Pouch Gen. 2 attaches not only to standard MOLLE gear but also to a wider range of tactical setups. Whether securing to belts, vests, or non-traditional gear, this pouch adapts to fit your needs while maintaining fast accessibility.

Learn more about this and other Products with a Mission® at www.NARescue.com



TECHNICAL SPECIFICATIONS

- Material**
Constructed with military-grade nylon and Cordura 1000, ensuring durability and resistance to harsh conditions.
- Reinforcement**
Features ipcabond Bar-Tack stitching on key stress points to ensure reliability during intense use.
- Dimensions**
Compact and lightweight, with customizable MOLLE straps and Velcro straps on the top and bottom, allowing easy attachment to a wide range of gear.



Tourniquets Sold Separately

DIMENSIONS

- Dimensions: L 5 in. x W 2 in. x D 1 in.
- Weight: 2.3 oz.

EZ-TOURNIQUET POUCH GEN2

ITEM #	COLOR
30-0559	Black
30-0560	Coyote
30-0561	Rgr Green

* Explore other Tourniquet holder options at www.NARescue.com

- Compatibility**
Specifically crafted to fit C-A-T® and similar tourniquet models for universal application.
- Attachment**
Single MOLLE, along with two sets of Velcro straps on the top and bottom, allowing easy attachment to a wide range of gear, including vests, belts, packs and more.
- Finishing**
Laser-cut cordura 1000, heat-sealed edges and 20mm seam allowance provide added strength and prevent fraying.

EZ-TOURNIQUET POUCH (GEN 2)



Tourniquet sold separately

