

USCG JOINT FIRST AID KIT (JFAK)



Learn more about this and other Products with a Mission® at www.NARescue.com

The Joint First Aid Kit (JFAK) from North American Rescue® is a low-profile, compact tactical medical kit designed for rapid response to life-threatening injuries. Built to meet the stringent requirements of the United States Coast Guard, this mission-ready JFAK integrates easily into any loadout—stow it inside an aid bag or mount it to body armor, assault vests, or packs using MOLLE/PALS-compatible webbing.

Configured to address the leading causes of preventable death in combat, the JFAK includes a C-A-T® Tourniquet, trauma dressings, Combat Gauze hemostatic, and gauze for massive hemorrhage control; a nasopharyngeal airway for airway management; and the Enhanced ARS for treatment of tension pneumothorax. Additional components—nitrile gloves, eye shield, trauma shears, and a combat casualty card—support efficient casualty care in austere and high-threat environments.

Durable, modular, and streamlined, this tactical JFAK ensures critical medical supplies are accessible when seconds matter—bridging the gap between point-of-injury care and advanced medical support.

FEATURES/SPECIFICATIONS


- Durable 500D Nylon
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Sleek Design with no Front Pocket

USCG JOINT FIRST AID KIT (JFAK)

ITEM # COLOR

80-1977	Black
---------	-------

KIT INCLUDES:



DIMENSIONS

- Case Closed: L 7 in. x W 7 in. x D 2.5 in.
- Weight: 1 lb. 5.8 oz.



KIT CONTENTS

- C-A-T® (Combat Application Tourniquet®) - Black
- QuikClot® Z-Fold Combat Gauze - 3 in. x 4 yd.
- NAR® S-Rolled Gauze - 4.5 in. x 4.1 yd.
- Flat Emergency Trauma Dressing (ETD®) - 4 in.
- Nasopharyngeal Airway with Lube - 28F
- HyFin® Compact (Vented) Dressing, Twin Pack
- GeckoGrip® Rescue Tape (6 pk)
- 2pr Nitrile Gloves - Large
- Black Permanent Marker - Large
- 1pr Mini Trauma XShears - Black/Black
- Enhanced ARS® Decompression Needle - 10ga x 3.25 in.
- PES Eye Shield
- Combat Casualty Card
- USCG JFAK bag - Black

