

# MARITIME GUARDIAN IFAK (BAG ONLY)



Learn more about this and other *Products with a Mission*® at [www.NARescue.com](http://www.NARescue.com)

## Rapid Access and Organization

Designed to solve problems, not create them, the trim USCG IFAK Bag (available in two of North American Rescue's most popular colors) can be tucked into small spaces or carried about with ease until it is required for service.

Available in NAR's 11 point Berry Amendment Nylon Standard.

## FEATURES/SPECIFICATIONS

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-TEX IR signature reduction hardware
- MOLLE-style Attachment System Compatible
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Sleek Design with no Front Pocket

## MARITIME GUARDIAN IFAK

ITEM #	COLOR
80-0150	Black
80-0157	COY

## DIMENSIONS

- H 8.25 in. X W 7.1 in X D 5.5 in.

