

# DELTA APRU ACCESSORIES

## Delta APRU Custom Foam Transport Case

The Delta Autonomous Portable Refrigeration Unit (APRU) Custom Foam Transport Case is a heavy duty plastic transport box with laser-cut, custom foam for optimal protection of the Autonomous Portable Refrigeration Unit (APRU) and its accessories.

### FEATURES/SPECIFICATIONS

- Heavy duty, with laser-cut custom foam inserts
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) and its accessories
- APRU and its accessories are sold separately

### DIMENSIONS

- Dimensions: H 30 in. x W 18 in. x D 23 in.
- Weight: 41 lbs.




*Contents Sold Separately*

### DELTA TRANSPORT CASE

ITEM #	NSN
35-0035	----

## DELTA APRU VEHICLE MOUNT & HARDWARE CABLE

ITEM #	NSN
35-0039	----



## Delta APRU Vehicle Mount & Hardware Cable

The Delta APRU Vehicle Mount and Hardwire Cable stabilizes the Autonomous Portable Refrigeration Unit (APRU) when transported in a vehicle. Includes hardwire cable and installation support. Does not include strap.

### FEATURES/SPECIFICATIONS

- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately
- Does not include strap

### DIMENSIONS

- Vehicle Mount Dimensions: H 1.25 in. x W 10 in. x L 19.25 in.
- Hardwire Cable Length: 15 ft.
- Total Weight: 5.2 lbs.



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# DELTA APRU ACCESSORIES

## Delta APRU Backpack Carrier with Pouches

The Delta APRU Backpack Carrier with Pouches is the CRO Medical (CRO) carrier for the Autonomous Portable Refrigeration Unit (APRU), adding carrying capability by hand or as a backpack. Also provides additional impact protection. Includes 2 extra battery pouches and a power cable storage pouch.

### FEATURES/SPECIFICATIONS

- Adds carrying capability by hand or as a backpack
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately
- Includes two extra battery pouches and a power cable storage pouch

### DIMENSIONS

- Dimensions: H 19 in. x W 9.5 in. x L 10.5 in.
- Weight: 3.2 lbs.



### DELTA BACKPACK CARRIER

ITEM #	NSN
35-0037	4110-01-706-0999

## DELTA SOFT SHOULDER CARRIER

ITEM #	NSN
35-0036	----



## Delta APRU Soft-Shoulder Carrier (FTECH)

The Delta APRU Soft-Shoulder Carrier is compatible with the Autonomous Portable Refrigeration Unit (APRU), adding carrying capability by hand. Also provides additional impact protection.

### FEATURES/SPECIFICATIONS

- Adds carrying capability by hand
- Provides additional impact protection
- Compatible with the Autonomous Portable Refrigeration Unit (APRU) - sold separately

### DIMENSIONS

- Dimensions: H 19.5 in. x W 10 in. x L 11 in.
- Weight: 3 lbs.



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

# DELTA APRU ACCESSORIES

## Delta APRU Battery-BB2590

The Delta APRU Battery - BB-2590 is a 7.5 Ah, rechargeable Lithium-Ion Battery. It is a replacement battery for the Autonomous Portable Refrigeration Unit (APRU), which is sold separately.

### FEATURES/SPECIFICATIONS

- Rechargeable Lithium-Ion Battery
- Compatible with the Delta Autonomous Portable Refrigeration Unit (APRU) - sold separately
- Designed and manufactured in the USA
- Nominal Voltage: 14.4V, 28.8V
- Ah: In 24V Mode: 7.5 Ah, In 12V Mode: 15
- Maximum Voltage: 16.8V
- Operating Temperature: -20°C to +60°C (-4°F to +140°F)
- Storage Temperature: -40°C to +40°C (-40°F to +104°F)
- Connector: Floating type per U.S. Army DWG # SC-C-179495
- State of Charge Indicator: Yes
- US Gov Spec: MIL-PRF-32052/1
- Battery Type: Rechargeable
- Smart Battery: SMBus

### DELTA APRU BATTERY - BB2590

ITEM #	NSN
35-0039	----



### DIMENSIONS

- Dimensions: H 5 in. x W 2.4 in. x L 4.4 in.
- Weight: 3.1 lbs.



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277