

NAR® K-9 LITTER

The NAR® K-9 Litter is a lightweight, compact litter designed for the movement of injured canines. The nylon litter can be packed to store in a cargo pocket or pouch and can be quickly deployed. The litter has leg cut outs to place the canine in the litter for safe movement. The K-9 Litter also has six large and four small carry handles for balanced carrying and allows for one to four person carries. The maximum weight capacity for each K-9 litter is 700 lbs.

**Learn more about this and
other *Products with a Mission*®
at www.NARescue.com**



FEATURES/SPECIFICATIONS

- Lightweight, compact litter for the movement of an injured canine
- Leg cut outs to safely carry the canine in the litter
- Six large and four small carry handles that allow for one to four person carries
- Can be carried in cargo pocket or pouch

K-9 LITTER

ITEM #

60-0074W

DIMENSIONS

- Litter: W 35 in. x H 48 in.
- Packaged: W 6 in. x H 8 in. x D 1 in.
- Weight: 15.3 oz
- Weight Capacity: 700 lbs



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277