

# CAT<sup>®</sup> HOLDER

## *Fits Securely and Easily on the Belt...*

The C-A-T<sup>®</sup> Holder was specifically designed to allow personnel to place their C-A-T<sup>®</sup> Tourniquet on their vest or gear for rapid and easy access.

It is made of 1000D IR Signature Reduced Nylon with MOLLE/PALS-style connectors.

The C-A-T<sup>®</sup> Holder protects your tourniquet from the elements and has an easy-open elastic pull tab. *Accommodates both Gen 6 and Gen 7 C-A-T<sup>®</sup> Tourniquets (not included)*

### Features:

- Tailored specifically for your C-A-T<sup>®</sup>, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- 1000D IR Signature Reduced Nylon
- Easy-open pull elastic tab with tourniquet ID patch – “TQ”
- NATO stock numbers available for COY and ODG. Contact Customer Service for details.

### Dimensions:

- L 6.25 in. x W 1.75 in. x D 1.5 in.
- Max belt width for horizontal mount: 1.5 in.
- Weight: 1.26 oz




REV010919



**NORTH AMERICAN RESCUE<sup>®</sup>**

www.NARescue.com • 888.689.6277

### Nylon Colors:



### C-A-T<sup>®</sup> Holder

\*NATO NSN

Color	Item #	NSN#
BLK	30-0057	
COY	30-0044	
RED	30-0128	
ODG	30-0053	6530-01-605-9432*
MTC	30-0059	
RTC	30-0207	
RGR	30-0208	