

NAVAL FIRST AID BOX RESPONSE KIT

Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels

Naval First Aid Box Response Kit Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels. Provides the first responder with rapid access and organization of critical life saving equipment for rendering on scene emergency care at sea. Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self-repairing zippers with silent 550 cord pulls. Lightweight, durable slimline pouch is securable using anti-pilferage seals & features a clear window to readily view contents . Pouch easily fits into the standard Naval First Aid Box and is ideal for restocking the Naval Gun Response Kit as needed.



Learn more about this and other **Products with a Mission®** at www.NARescue.com



Pouches easily fit into the standard Naval First Aid Box (not included)

FEATURES/SPECIFICATIONS

- Designed to meet the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels
- Constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self-repairing zippers with silent 550 cord pulls.
- Lightweight, durable slimline pouch is securable using anti-pilferage seals & features a clear window to readily view contents
- Pouch easily fits into the standard Naval First Aid Box and is ideal for restocking the Naval Gun Response Kit as needed.

NAVAL FIRST AID BOX RESPONSE KIT

| ITEM # | OPTION | NSN |
|---------|-----------|--|
| 80-0419 | Kit | 6545-01-648-2705 <small>R_X ONLY</small> |
| 80-0357 | Bags Only | 8465-01-666-7328 |

DIMENSIONS

- L 15 in. x H 6 in. x D 4 in.
- Weight: 1 lb 15 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

NAVAL FIRST AID BOX RESPONSE KIT

all KITS INCLUDE:



KIT CONTENTS

- 1x Pouch (CCRK® - Naval First Aid Box Response)
- 1x C-A-T® (Combat Application Tourniquet®)
- 1x ETD® 6 in. Emergency Trauma Dressing
- 1x ETD® Abdominal Emergency Trauma Dressing
- 1x Nasopharyngeal Airway 28F
- 1x Nasopharyngeal Airway 34F
- 1x Surgical Lubrication
- 1x CPR Mask, Microshield
- 1x HyFin® Vent Chest Seal Twin Pack
- 1x Fox Eye Shield with Garter
- 1x Roll 2 in. Surgical Tape
- 2x SAM® Splint II
- 1x Elastic Wrap Bandage, 3 in.
- 1x Trauma Shears (7.25 in.)
- 1 pr Bear Claw® Nitrile Gloves, Large
- 1x Permanent Marker, Large



NAVAL FIRST AID BOX RESPONSE KIT

| ITEM # | OPTION | NSN |
|---------|-----------|-------------------------------------|
| 80-0419 | Kit | 6545-01-648-2705 ^{RX} ONLY |
| 80-0357 | Bags Only | 8465-01-666-7328 |

Learn more about this and other **Products with a Mission®** at www.NARescue.com

