

PATROL VEHICLE TRAUMA BAG ONLY

Learn more about this and other Products with a Mission® at www.NARescue.com



The Patrol Vehicle Trauma Bag from North American Rescue is a modular storage bag designed to house essential medical equipment in an easily accessible, quick to deploy platform ideal for carriage in your vehicle. This durable zippered nylon bag features two separate interior compartments providing optimal organization of your gear in any number of custom configurations. The bag's exterior also houses additional equipment storage such as a convenient sheath for your trauma shears and an outer storage pocket on one side and a MOLLE-style system on the opposite side to accommodate extra equipment as needed, such as additional tourniquets, ammo, or personal protective equipment. This is the bag only. All equipment is sold separately.



DIMENSIONS

- Closed: H 7 in. x W 12 in. x D 5 in.
- Weight: 12.8 oz

FEATURES/SPECIFICATIONS

- Rugged nylon zippered bag with easy to carry handles
- Easily accessible, quick to deploy platform for a medical kit in any vehicle
- Modular storage system features 2 separate interior compartments
- Bag only. All equipment is sold separately

for optimal organization of gear

- Exterior includes MOLLE-style system on one side and trauma shears sheath/extra storage pocket on the other to accommodate easy to access equipment



PATROL VEHICLE TRAUMA BAG ONLY

ITEM #	COLOR	NSN
80-0531T	Black	----
80-1970T	Red	----

