

NAR Simple Thoracostomy Kit

The NAR Simple Thoracostomy Kit provides the essential equipment needed to perform a simple (finger) thoracostomy. Packaged in a compact and easy-to-open vacuum-sealed bag, the kit contains all necessary components to perform the procedure in one package. Contents include a 3-pack Chloraprep® swapstick for skin prep, a #10 PenBlade protected scalpel, sterile gauze, curved Rochester Pean hemostat to penetrate parietal pleura and dissect intercostal muscle, HyFin Compact Twin Pack Chest Seal for spontaneously breathing patients and a marker to write on chest seal as needed.

FEATURES/SPECIFICATIONS

- Compact kit that contains all necessary equipment to perform simple thoracostomy
- Includes 3-pack skin prep for performing the procedure bilaterally if needed
- HyFin Compact Vented Twin Pack Chest Seal to cover site if required

EDUCATIONAL TRAINING RESOURCES AVAILABLE

Download **NAR Simple Thoracostomy Kit Training Power Point Slides**



View online **NAR Simple Thoracostomy Kit Product Training Video**




Learn more about this and other **Products with a Mission®** at www.NARescue.com



DIMENSIONS

- Packaged: L 4.5 in. x W 7.5 in. x D 0.75 in.
- Weight: 5.28 oz

KIT CONTENTS

- 1x Rochester Pean Curved - 6.25 in.
- 1x Protected Scalpel - #10
- 1x HyFin® Vent Compact Chest Seal -Twin Pack
- 2x Gauze Pads (2pk) - 4 in. x 4 in.
- 1x Chloraprep Swabsticks (3pk)
- 1x Permanent Marker - Small

SIMPLE THORACOSTOMY KIT

ITEM #	NSN #
85-3120	6515-01-733-0892

R_X ONLY

SIMPLE THORACOSTOMY KIT CONTENTS:

QTY 1 Hemostat 6.25" Rochester Pean Curved

QTY 1 Protected #10 Scalpel

QTY 1 HyFin® Vent Compact,Twin

QTY 2 Gauze Pads 4x4"(2pk)

QTY 1 Chloraprep Swabsticks,3pk

QTY 1 Permanent Marker, Small