

NAR-4™ AID KIT

The constant challenge of striking a balance between the increasing amount of equipment carried into battle and the need to optimize performance is the impetus behind this low-profile, load carriage platform for Tactical Healthcare Professionals.

The NAR-4™ Aid Kit is an ergo-dynamic, durable, functional solution that supports Tactical Combat Casualty Care equipment standards. The NAR-4™ features a modular system with 5 removable nylon bags that can be configured to carry airway tools, IV equipment & hemorrhage control gear.

The slimline outer bag can be carried as a low-profile backpack using shoulder straps which are stowed away until needed or pre-rigged to attach directly to body armor using four rapid attachment straps.



Learn more about this and other Products with a Mission® at www.NARescue.com



Ergo-Dynamic, Durable and Functional Solution for Tactical Healthcare Professionals

FEATURES/SPECIFICATIONS

- Configurable based on individual skill set and mission requirements
- Rapid access and organization is enhanced by segregated removable equipment pouches
- 500D CORDURA® Berry Amendment Nylon Bag with GhillieTEX™ hardware
- Sand and water resistant, self-repairing zippers with silent 550 cord pulls
- PALS/MOLLE-style attachment system on front and side
- Packing list is consistent with TCCC Guidelines

for a reduced IR signature



NAR-4 AID KITS

ITEM #	COLOR	HEMOSTATIC	NSN
80-0309	Black	ChitoGauze	----
80-0310	COY	ChitoGauze	----
80-0313	MTC	ChitoGauze	----
80-0320	ODG	ChitoGauze	----
80-0181	Black	Combat Gauze	----
80-0183	COY	Combat Gauze	----
80-0244	MTC	Combat Gauze	----
80-0308	ODG	Combat Gauze	----


DIMENSIONS

- Cased Closed:
H 19 in. X W 15 in. X D 7 in.
- Weight: 12 lb. 8.3 oz.



NAR-4™ AID KIT

all KITS INCLUDE:



KIT CONTENTS

- 2- C-A-T™ (Combat Application Tourniquet®)
- 2- NAR® Emergency Trauma Dressing, 6 in.
- 1- NAR® Abdominal Emergency Trauma Dressing
- 2- Choice of Hemostatic (Combat Gauze® or ChitoGauze® Pro)
- 4- NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1- King LTS-D™ Supraglottic Airway Kit
- 2- NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1- Tactical Suction Device
- 1- HyFin® Vent Chest Seal Twin Pac
- 2- ARS® for Needle Decompression (14G x 3.25 in.)
- 1- Cyclone® New & Improved BVM
- 1- BOA® IV Constricting Band
- 4- NAR® Saline Lock Ki
- 1- NAR® Sharps Shuttle
- 1- Armadillo™ Medication Storage Case
- 1- NAR® Trauma Shears (7.25 in.)
- 2- GeckoGrip™ Multipurpose Tape
- 2- SAM® Splint II
- 1- CT-6 Traction Splint
- 1- Tactical CrickKit®
- 6- Polycarbonate Eye Shield - Orange
- 1- Permanent Marker - Black, Lg
- 8 pr- Black Talon® Nitrile Trauma Gloves
- 4- Combat Casualty Card (TCCC Triage Card)
- 1- Combat Casualty Reference Card
- 1- Bag (CCRK® - NAR-4™)

Learn more about this and other Products with a Mission® at www.NARescue.com



NAR-4 AID KITS

ITEM #	COLOR	HEMOSTATIC	NSN
80-0309	Black	ChitoGauze	----
80-0310	COY	ChitoGauze	----
80-0313	MTC	ChitoGauze	----
80-0320	ODG	ChitoGauze	----
80-0181	Black	Combat Gauze	----
80-0183	COY	Combat Gauze	----
80-0244	MTC	Combat Gauze	----
80-0308	ODG	Combat Gauze	----

