

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

# HEAT REFLECTIVE SHELL (HRS)



## Excellent Thermal Insulation...

The reinforced Heat Reflective Shell (HRS™) from North American Rescue® is a strong, flexible, lightweight covering designed to completely encapsulate the patient during casualty care. The HRS™ is constructed of a 4-ply, composite fabric with a protected non-conductive reflective layer that provides excellent thermal insulation and durable support for up to 250+ lb. Rapid 360° access to the casualty with minimal exposure to the elements is achieved through the use of 1.5 in. continuous Velcro® closures. The HRS™ has a built-in hood, integrated fluid absorption pad and top-to-bottom tapered shape to maximize isothermal capabilities.

## DIMENSIONS

- Packaged: H 10 in. x W 6.75 in. x D 2.5 in.
- Size: L 78 in. x W 43 in.
- Weight: 1 lb 13 oz



### Heat Reflective Shell

[REF] ZZ-0051 [REF] ZZ-0051-EU NSN: 6532-01-627-7888

U.S. Patent 7,766,950 B2

USA ONLY  
NOT MADE WITH NATURAL RUBBER LATEX

EC REP CE

Assembled in Mexico from U.S. and foreign components  
REV090221

**NORTH AMERICAN RESCUE**  
www.NARescue.com • 888.689.6277  
35 Tedwall Court • Greer, SC 29650 • USA

LOT

Intended Use: The device is intended to reflect heat and provide thermal insulation.

In the event of a malfunction follow local protocols and report any serious incident to North American Rescue or authorized representative and the competent authority for the member state.

Within the European Union (EU), the device is not indicated for pregnant and pregnant women due to EU regulatory requirements.

MD

**1.**

**2.**

**3.**

## FEATURES/SPECIFICATIONS

- » 4-ply composite fabric with protected non-conductive reflective layer that provides excellent thermal insulation
- » Tapered shape from top to bottom to maximize isothermal capabilities
- » Built-in hood & integrated fluid absorption pad

- » Continuous 1.5 in. Velcro® closures for rapid 360° access to the casualty with minimal exposure to the elements
- » Reinforced and weight tested to support 250+ lb
- » Vacuum-packed for small cube footprint
- » U.S. Patent 7,766,950 B2

## HEAT REFLECTIVE SHELL (HRSI)

ITEM #	NSN
ZZ-0051	6532-01-627-7888
ZZ-0051-EU	6532-01-627-7888

**NORTH AMERICAN RESCUE**  
www.NARescue.com • 888.689.6277