

NASOPHARYNGEAL AIRWAY WITH LUBRICANT

The Nasopharyngeal Airway 26F with Lubricant shows North American Rescue's attention to details. This NPA was designed to be dependable, efficient, comfortable, and secure.

Efficient. The product includes a conveniently attached sealed packet of surgical lubricant to facilitate easy insertion.

Designed to be more compatible with human anatomy than the standard oropharyngeal airway, this product is recommended for use as a nasal airway adjunct in semi-conscious or unconscious patients with an intact gag reflex.

Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Designed for semi-conscious or unconscious patients with an intact gag reflex
- Packaged with surgical lubricant attached

DIMENSIONS

- Packaged: L 6 in. x W 16 in. x D 0.5 in.
- Weight: 0.5 oz



NASOPHARYNGEAL AIRWAY W/ LUBE

ITEM #	SIZE	LENGTH	NSN
ZZ-1360	20F	105 mm	----
ZZ-0424	22F	105 mm	----
ZZ-0422	24F	105 mm	6515-01-687-1224
ZZ-0387	26F	110 mm	----
ZZ-0034	28F	117 mm	6515-01-529-1187
ZZ-0423	30F	118 mm	----
ZZ-0388	32F	140 mm	6515-01-687-1244
ZZ-0776	34F	155 mm	----
ZZ-0777	36F	170 mm	----

* There are other options available at a place for this item and other items.

RX ONLY



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277