

# NAR-4™ CHEST POUCH

The auxiliary NAR-4™ companion is streamlined to provide Tactical Healthcare Professionals with rapid access to critical basic & advanced equipment on their chest.

The container is constructed of 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction, GhillieTEX™ IR signature reduction hardware and a silent, self repairing nylon zipper with 550 pull cords.

The kit is designed to be secured to the front of the body armor using a MOLLE-style attachment system.

The NAR-4™ Chest Pouch works seamlessly with the NAR® Quick Detach Panel (sold separately)



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



*Streamlined to provide Tactical Healthcare Professionals rapid access...*

## FEATURES/SPECIFICATIONS

- Configurable based on individual skill set and mission requirements
- Rapid access and organization is enhanced by elastic storage loops that keep items in place until needed
- 500D CORDURA® Solution-Dyed Nylon with built in Near Infrared Signature Reduction
- Sand and water resistant zippers with silent pulls
- PALS/MOLLE-style attachment system
- Ambidextrous Trauma Shear configuration
- Compatible with NAR® Quick Detach Panels (sold separately)



## NAR-4™ CHEST POUCH

ITEM #	COLOR	HEMOSTATIC	NSN
80-0283	Black	ChitoGauze	----
80-0284	COY	ChitoGauze	----
80-0287	MTC	ChitoGauze	----
80-0288	ODG	ChitoGauze	----
80-1664	Red	ChitoGauze	----
80-0173	Black	Combat Gauze	----
80-0175	COY	Combat Gauze	----
80-0247	MTC	Combat Gauze	----
80-0289	ODG	Combat Gauze	----
80-1665	ODG	Combat Gauze	----

## DIMENSIONS

- Cased Closed: L 7.75 in. x W 9.75 in. x 4.5 in.
- Weight: 2 lb. 10.3 oz.



# NAR-4™ CHEST POUCH

**all** KITS INCLUDE:



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## KIT CONTENTS

- 2- C-A-T® (Combat Application Tourniquet®) (14G x 3.25 in.)
- 2- NAR® Emergency Trauma Dressing, 6 in.
- 2- NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 2-Choice of Hemostatic (Combat Gauze® or ChitoGauze® Pro)
- 2- NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1- Tactical CricKit®
- 1- HyFin® Vent Chest Seal Twin Pac
- 2- ARS® for Needle Decompression
- 3- Black Talon® Nitrile Glove Kit, Lg.
- 2- Combat Casualty Card
- 1- Scissor Leash™
- 1- NAR® Trauma Shears (7.25 in.)
- 1- Permanent Marker
- 2- Polycarbonate Eye Shield
- 1- Bag (CCRK® - NAR-4™ Chest Pouch)



## NAR-4™ CHEST POUCH

ITEM #	COLOR	HEMOSTATIC	NSN
80-0283	Black	ChitoGauze	----
80-0284	COY	ChitoGauze	----
80-0287	MTC	ChitoGauze	----
80-0288	ODG	ChitoGauze	----
80-1664	Red	ChitoGauze	----
80-0173	Black	Combat Gauze	----
80-0175	COY	Combat Gauze	----
80-0247	MTC	Combat Gauze	----
80-0289	ODG	Combat Gauze	----
80-1665	ODG	Combat Gauze	----