

# TECHLINE TRAUMA WEARABLE THIGH LACERATION WOUND

## Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma Wearable Thigh Laceration simulates a deep laceration exposing fat and muscle tissue with tourniquet and pressure responsive blood flow.

The Techline Trauma® Wearable Thigh Laceration Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world deep laceration. The soft, flexible, light toned skin wraps comfortably around the thigh and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue and muscle and includes an integrated bleeding system tube that allows role players to demonstrate active bleeding during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Thigh Laceration is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.

Role player use requires blood supply pump (sold separately).



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Tourniquet/Direct Pressure Responsive
- Integrated bleeding system
- Durable and UV Resistant

## WEARABLE THIGH LACERATION WOUND

ITEM #

93-0203

## DIMENSIONS

- Dimensions: L 10 in. x W 9 in. x H 1.5 in.
- Weight: 2 lbs.

\* Explore other Techline Trauma products at [www.NARescue.com](http://www.NARescue.com)