

# TECHLINE TRAUMA WEARABLE SEVERE LEFT ARM BURN

## Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma® Wearable Severe Left Arm Burn simulates third degree burns to the left arm. Wrapping comfortably around the actor or manikin's left arm simulated and attached with hook and loop fastening.

The Techline Trauma® Wearable Severe Left Arm Burn Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world third degree burn injury. The soft, flexible, light toned skin wraps comfortably around the left forearm and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue and muscle that allows role players or manikins to demonstrate third degree burns during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Severe Burn Arm Left is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.



Learn more about this and other **Products with a Mission®** at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Durable and UV Resistant
- Approximate burn of ~4% of TBSA

## WEARABLE SEVERE LEFT ARM BURN

ITEM #

93-0224

\* Explore other Techline Trauma products at [www.NARescue.com](http://www.NARescue.com)

## DIMENSIONS

- Dimensions: L 20 in. x W 6 in. x H 3 in.
- Weight: 2 lbs.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277