

# TECHLINE TRAUMA WEARABLE PACKABLE EVISCERATION

## Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma® Wearable Packable Evisceration simulated an evisceration exposing the bowels and tissue, offering the opportunity to treat a traumatic abdominal protrusion.

The Techline Trauma® Wearable Packable Evisceration Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world evisceration injury. The soft, flexible, light toned skin wraps comfortably around the abdomen and attaches with the sewn hook and loop connectors. This realistic injury features exposed bowel and fat tissue and includes an integrated bleeding system tube that allows role players or manikins to demonstrate active bleeding during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Evisceration is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.

Role player use requires blood supply pump (sold separately).



Learn more about this and  
other **Products with a Mission®**  
at [www.NARescue.com](http://www.NARescue.com)



## FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Integrated bleeding system
- Packable Bowel can be fully inserted, or protrude depending on training requirements
- Durable and UV Resistant

## WEARABLE PKBL EVISCERATION

ITEM #

93-0080

\* Explore other Techline Trauma  
products at [www.NARescue.com](http://www.NARescue.com)

## DIMENSIONS

- Dimensions: L 15 in. x W 14 in. x H 1 in.
- Weight: 5 lbs.



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277