

NAR® NASOPHARYNGEAL AIRWAY

28F (117 MM) - NON-LUBRICATED

Introducing the North American Rescue® Non-Lubricated Nasopharyngeal Airway




Our non-lubricated nasopharyngeal airway is an airway adjunct that can be used on semi-conscious or unconscious patients with or without an intact gag reflex. The beveled tip and non-latex, soft flexible material avoids damage to the nasal passage. Easy to use and packaged in a durable pouch with low cube size and weight. The non-lubricated 28F (117mm) nasal airway is anatomically shaped and made of soft material to offer optimal conformity, reducing the risk of tissue damage and enhancing patient safety. **Made in USA** with high grade materials and designed for optimal efficiency and ease of use.

Learn more about this and other **Products with a Mission®** at www.NARescue.com



FEATURES/SPECIFICATIONS

- Single Nasopharyngeal device
- Rounded bevel tip
- Soft material avoids damage to nasal passage
- Unique proprietary material and design maintains patent shape
- Durable packaging
- Low cube and weight
- Sterile
- Latex-free
- Made in USA



NAR® NASOPHARYNGEAL AIRWAY, 28F - NON-LUBRICATED

QTY	ITEM #	NSN
Single	10-0061	----
Box of 10	10-0068	6515-01-710-9848

DIMENSIONS

- Ind. Package: L 7.5 in. x H 3.5 in.
- Box: L 9 in. x W 4 in. x D 4 in.
- Device Size: 28 FR (117 mm)
- Weight: 0.02 lbs (Ind.)
0.50 lbs (Box)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277