

TECHLINE TRAUMA WEARABLE MULTIPLE GSW ABDOMEN

Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma® Wearable Multiple GSW Abdomen Simulator represents a through and through gunshot wound exposing fat and muscle. All exit wounds are packable, one entry and exit can be packed through and through. This wearable wound bleeds from all entry and exits with tourniquet and pressure responsive blood flow.

The Techline Trauma® Wearable Multiple GSW Abdomen Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world abdominal GSW injury. The soft, flexible, light toned skin wraps comfortably around the abdomen and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue, fat and muscle and includes an integrated bleeding system tube that allows role players or manikins to demonstrate active bleeding during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Multiple GSW Abdomen is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.

Role player use requires blood supply pump (sold separately).

Learn more about this and other Products with a Mission® at www.NARescue.com



DIMENSIONS

- Dimensions: L 18 in. x W 14 in. x H 1 in.
- Weight: 4 lbs.

FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Packable Sites
- Direct Pressure Responsive
- Integrated bleeding system
- Durable and UV Resistant

MULTIPLE GSW ABDOMEN

ITEM #

93-0229

* Explore other Techline Trauma products at www.NARescue.com



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277