

# TECHLINE TRAUMA WEARABLE IMPALED OBJECT WOUND

## Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma® Wearable Impaled Object represents a piece of rebar impaled into the flesh. With a range of applicable areas on a role player or manikin, this simulator offers variety to any training program.

The Techline Trauma® Wearable Impaled Object Simulator from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device enabling live action role players to simulate a real-world impalement injury. The soft, flexible, light toned skin wraps comfortably around the arm or leg and attaches with the sewn hook and loop connectors. This realistic injury features exposed tissue and muscle that allows role players or manikins to demonstrate impalement during scenarios. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Impaled Object is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.



**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



## DIMENSIONS

- Dimensions: L 8 in. x W 8 in. x H 4 in.
- Weight: 1 lb.

## FEATURES/SPECIFICATIONS

- Hook and loop closure
- Exposed tissue and muscle
- Durable and UV Resistant
- Multiple Site Application

## IMPALED OBJECT WOUND

ITEM #

93-0205

\* Explore other Techline Trauma products at [www.NARescue.com](http://www.NARescue.com)



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277