

# TECHLINE TRAUMA WOUND EXTENDER

## Minimizes Moulage Requirements – While Maximizing Realism.

The Techline Trauma® Wearable Wound Extender from North American Rescue® Simulation and Techline Technologies, Inc., is a silicone prosthetic device allowing placement of Techline Wearable wounds around live action role players and manikins of different shapes and sizes. The soft, flexible, light toned skin wraps comfortably around a variety of locations and attaches with the sewn hook and loop connectors. Enhancing realism creates a valuable training event that better prepares first responder trainees for real-world situations. Easily cleaned and reusable, this tool enables rapid preparation for training exercises and minimizes moulage requirements – while maximizing realism.

The Techline Trauma® Wearable Wound Extender is manufactured by Techline Trauma® and is compatible with all TOMManikin® and TAMlkin® simulators.



**Impalement Wound sold separately**

**Learn more about this and  
other Products with a Mission®  
at [www.NARescue.com](http://www.NARescue.com)**



## DIMENSIONS

- Dimensions: L 9 in. x W 7 in. x H 1 in.
- Weight: 1 lb.



## FEATURES/SPECIFICATIONS

- Hook and loop closure
- Extends Wearable Wounds for Placement on a Variety of Role Player or Manikin Locations or Sizes
- Promotes Versatility of Wound Placement
- Durable and UV Resistant
- Wound Extender ONLY. Other wounds sold separately

## WOUND EXTENDER

ITEM #

93-0214

*Explore other Techline Trauma products at [www.NARescue.com](http://www.NARescue.com)*



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277