

NAR® SIMPLE THORACOSTOMY 2.0 KIT

Learn more about this and other Products with a Mission® at www.NARescue.com



The NAR® Simple Thoracostomy Kit - 2.0 provides the essential equipment needed to perform a simple (finger) thoracostomy. Packaged in a compact and easy-to-open vacuum-sealed bag, the kit contains all necessary components to perform the procedure in one package. Contents include two 3-pack Chloraprep® swabsticks for skin prep, a #10 protected safety scalpel, sterile gauze, curved Rochester Pean hemostat to penetrate parietal pleura and dissect intercostal muscle, a sterile fenestrated drape to create a sterile barrier, and a permanent marker to use as needed.

SIMPLE THORACOSTOMY 2.0 KIT CONTENTS:



KIT CONTENTS

- 1 x #10 Protected Scalpel
- 1 x Curved Rochester Pean Hemostat 6.25 in.
- 1 x Sterile Fenestrated Drape - 18 in. x 25 in.
- 2 x Gauze Pad (2pk) - 4 in. x 4in.
- 2 x Chloraprep® Swabsticks (3pk)
- 1 x Surgical Tape - 1 in.
- 1 x Black Permanent Marker - Small

DIMENSIONS

- Packaged: L 6 in. x W 9 in. x D 1.5 in.
- Weight: 0.39 lbs. (6.17 oz.)

FEATURES/SPECIFICATIONS

- Compact kit that contains all necessary equipment to perform simple thoracostomy
- Includes two (3-pack) skin prep swabsticks for performing the procedure bilaterally if needed

NAR SIMPLE THORACOSTOMY 2.0 KIT

ITEM #	NSN
83-0040	6515-01-738-0746

Rx ONLY