

M-FAK

MINI FIRST AID KIT

Our most compact, half size First Aid Kit empowering First Responders with the most requested point-of-wounding equipment in the smallest cube space possible

The M-FAK, North American Rescue's Mini First Aid Kit, was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K. platform, the M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty. Treatment pack out options include Basic, Basic with Combat Gauze and Advanced versions.

The M-FAK's super compact, 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T.[®] tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips to attach additional gear such as a TQ/Holder (sold separately) if needed.



Shown: Advanced - ODG

DIMENSIONS

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz

FEATURES/SPECIFICATIONS

- NAR's most compact, versatile Mini First Aid Kit
- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Treatment pack out option: BLS (Basic Life Support)
- Super compact, rugged nylon platform that allows attachment both vertically and horizontally
- Half the size of our popular T.O.R.K. platform
- Compatible with NAR[®] Quick Detach Panels (sold separately)



Learn more about this and other Products with a Mission[®] at www.NARescue.com



M-FAK

ITEM #	COLOR	OPTION	NSN
80-0494	Black	Basic	----
80-0495	Coyote	Basic	----
80-0872	Gray	Basic	----
80-0497	OD Green	Basic	----
80-0996	RGR Green	Basic	----
80-0496	Multicam	Basic	----
80-0890	Red	Basic	----
80-0498	Black	with Combat Gauze [®]	----
80-0499	Coyote	with Combat Gauze [®]	----
80-0873	Gray	with Combat Gauze [®]	----
80-0501	OD Green	with Combat Gauze [®]	----
80-0998	RGR Green	with Combat Gauze [®]	----
80-0500	Multicam	with Combat Gauze [®]	----
80-0891	Red	with QuikClot [®]	----
80-0502	Black	Advanced**	----
80-0503	Coyote	Advanced*	----
80-0874	Gray	Advanced**	----
80-0505	OD Green	Advanced**	----
80-0997	RGR Green	Advanced**	----
80-0504	Multicam	Advanced**	----
80-0892	Red	Advanced**	----

* Explore other IFAK's at www.NARescue.com

** R_x ONLY



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

M-FAK

MINI FIRST AID KIT

all KITS INCLUDE:



QTY 1
C-A-T®
Tourniquet



QTY 1
ETD® 4 in.
Dressing, Flat





QTY 1
HyFin® Vent
Compact, Twin Pk



PAIR 1
Nitrile
Gloves, Lg

+plus, by kit type:



KIT CONTENTS

- 1x C-A-T® - Black (Combat Application Tourniquet®)
- 1x 4 in. Flat ETD®

- 1x HyFin® Vent Compact Chest Seal - Twin Pack
- 1 pr Bear Claw® Nitrile Trauma Gloves, lg.

+plus, by kit type:

- 1x S-rolled Gauze (4.5 in. x 4.1 yd) -OR-
- 1x Combat Gauze® LE Hemostatic Bandage -OR-
- 1x QuikClot® Bleeding Control Dressing* -OR-
- 1x Combat Gauze® LE Hemostatic Bandage
- 1x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 1x ARS® Needle Decompression Kit (14 G x 3.25 in.)

* Only available in Red.



Shown:
Basic - RGR

Learn more about this and other Products with a Mission® at www.NARescue.com

