

# RIG SERIES EAGLE IFAKS

## Low Profile IFAK with Rapid Access To Your C-A-T® Tourniquet

The RIG Series Eagle IFAK is a new low profile individual first aid kit designed with a protective hood and quick latch release to keep your life-saving equipment protected and rapidly accessible for deployment. Highlighted design advancements include easy-one hand quick release hood latch and zippers, detachable inner shock cord panel that detaches from the IFAK bag and attachment system that allows both vertical (MOLLE backing) and horizontal (3 in. belt loop) carriage orientations. Kits come in 3 configurations: Basic, with Combat Gauze, and Advanced.



Shown: Advanced - BLK

## FEATURES/SPECIFICATIONS

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Kits come in 3 configurations: Basic, with Combat Gauze, and Advanced
- Protective hood with quick latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Quick pull latch and tab for easy one-handed opening
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components & is easily removable for rapid access
- Compatible with NAR® Quick Detach Panels (sold separately)

## DIMENSIONS

- Packaged: H 8 in. x W 4 in. x D 3 in.
- Weight: 1.18 lb

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## M-FAK

ITEM #	COLOR	OPTION	NSN
80-0509	Black	Basic	----
80-0510	Coyote	Basic	----
80-0511	Multicam	Basic	----
80-0991	RGR Green	Basic	----
80-0517	Black	with Combat Gauze®	----
80-0518	Coyote	with Combat Gauze®	----
80-0519	Multicam	with Combat Gauze®	----
80-0993	RGR Green	with Combat Gauze®	----
80-0513	Black	Advanced**	----
80-0514	Coyote	Advanced**	----
80-0515	Multicam	Advanced**	----
80-0992	RGR Green	Advanced**	----

\* Explore other IFAK's at [www.NARescue.com](http://www.NARescue.com)

\*\* R<sub>X</sub> ONLY



# RIG SERIES EAGLE IFAKS

**all KITS INCLUDE:**



**QTY 1**  
C-A-T®  
Tourniquet



**QTY 1**  
ETD® 4 in.  
Dressing, Flat



**QTY 1**  
NAR® Gauze,  
Wound Packing




**QTY 1**  
HyFin® Vent  
Compact, Twin Pk



**PAIR 1**  
Black Talon  
Nitrile Gloves, LG

**+plus, by kit type:**



## KIT CONTENTS

- 1x C-A-T® - Black (Combat Application Tourniquet®)
- 1x NAR® Wound Packing Gauze
- 1x 4 in. Flat ETD®
- 1x HyFin® Vent Chest Seal Twin Pack
- 1 pr Black Talon® Nitrile Trauma Gloves, lg.

**+plus, by kit type:**

- 1x NAR® Wound Packing Gauze - OR -
- 1x Combat Gauze® LE Z-fold Hemostatic - OR -
- 1x Combat Gauze® LE Z-fold Hemostatic
- 1x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 1x ARS® Needle Decompression Kit (14 G x 3.25 in.)



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

