

# REFLEX IFAK SYSTEM KIT

Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)



## Carry organized medical supplies with ambidextrous accessibility

The Reflex™ IFAK System Kit is a two-piece system, designed to carry organized medical supplies with ambidextrous accessibility. This system, developed with direct input from active-duty medical personnel, is constructed primarily with heavy-duty nylon laminate which allows rapid deployment of medical supplies.

Reflex™ IFAK System Kits feature premium life saving medical components from North American Rescue that address the leading causes of preventable combat death identified by Tactical Combat Casualty Care (TCCC) protocols. These critical tools enable responders to stop major bleeding, address penetrating trauma, and to help prevent hypothermia when providing care at the point of wounding prior to transport to definitive care facilities. The system is available in three levels of care: Basic, Basic with Combat Gauze, and Advanced.



## REFLEX IFAK SYSTEM KIT

ITEM #	COLOR	OPTION	NSN
80-1127	Black	Basic	----
80-1128	Blue	Basic	----
80-1129	Coyote	Basic	----
80-1130	Gray	Basic	----
80-1131	MCB	Basic	----
80-1132	MTC	Basic	----
80-1133	OD Green	Basic	----
80-1134	Black	Basic with Combat Gauze®	----
80-1135	Blue	Basic with Combat Gauze®	----
80-1136	Coyote	Basic with Combat Gauze®	----
80-1137	Gray	Basic with Combat Gauze®	----
80-1138	MCB	Basic with Combat Gauze®	----
80-1139	MTC	Basic with Combat Gauze®	----
80-1140	OD Green	Basic with Combat Gauze®	----
80-1141	Black	Advanced	----
80-1142	Blue	Advanced	----
80-1143	Coyote	Advanced	----
80-1144	Gray	Advanced	----
80-1145	MCB	Advanced	----
80-1146	MTC	Advanced	----
80-1147	OD Green	Advanced	----

## FEATURES/SPECIFICATIONS

- Features HANK™ (High Abrasion-Resistant Neoprene Kevlar® composite) grab handles for long-term durability
- Can mount vertically or horizontally on belts 1.5 in.- 2.25 in., and horizontally on MOLLE
- Intentional organization features put all supplies within easy reach
- Carry handle and shock cord loops allow the roll to be attached to a patient, or hung from an I.V. rack
- "Mini" MOLLE allows for precise fit on belts as narrow as 1.5 in. or two-row MOLLE panels

## DIMENSIONS

- Measurements (English): L 8 in. x W 3 in. x H 5.5 in.
- Measurements (Metric): L 20.32 cm x W 7.62 cm x H 13.97 cm
- Weight (English): 1.58 lb
- Weight (Metric): 716.68 g

# REFLEX IFAK SYSTEM KIT

**all KITS INCLUDE:**




**+plus, by kit type:**



## KIT CONTENTS

- 1x C-A-T® - Black (Combat Application Tourniquet®)
  - 1x ETD®, 4 in. Flat (Emergency Trauma Dressing)
  - 1x HyFin® Vent Chest Seal Twin Pack
  - 1x NAR® Survival Blanket
  - 1x Trauma Shears (7.25 in.)
  - 1x NAR® Flat Tape
  - 2 pr Responder Black Nitrile Gloves, LG
  - 2x Combat Casualty Document Tool Card
  - 1x Permanent Marker - Large / Black
- +plus, by kit type:**
- 2x NAR® Wound Packing Gauze (3 in. x 5 yd.) - or -
  - 1x Combat Gauze Hemostatic Bandage (with Z-Fold)
  - 1x NAR® Wound Packing Gauze (3 in. x 5 yd.) - or -
  - 1x Combat Gauze Hemostatic Bandage (with Z-Fold)
  - 1x NAR® Wound Packing Gauze (3 in. x 5 yd.)
  - 1x NAR® Pre-Lubricated Nasopharyngeal Airway - 28F(117 mm)
  - 1x ARS® Needle Decompression Kit (10 G x 3.25 in.)



Learn more about this and other Products with a Mission® at [www.NARescue.com](http://www.NARescue.com)

