

EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

The Expeditionary Casualty Response Kit (ECRK) is a lightweight platform designed to rapidly access and organize all your essential trauma management equipment. Features include silent, dirt and sand-resistant zippers that close over a large main compartment with six (6) segregated removable equipment pouches; 4 Bleed Pouches, 1 Airway Pouch, and 1 Miscellaneous Pouch. Additionally, a separate lid compartment secures the R.E.E.L. Non-Rigid Litter while two internal pouches on the internal lid compartment and two external pouches on each side of the bag accommodate additional trauma management essentials.



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Lightweight platform to rapidly access and organize all your essential trauma management equipment
- Featuring a large main compartment with six segregated removable pouches secured by hook and loop
 - » 4 Internal Bleed Pouches with Red pull handle secured by hook and loop
 - » 1 Airway Pouch with Blue pull handle secured by hook and loop
 - » 1 Miscellaneous Pouch with Coyote pull handle secured by hook and loop
- Constructed of 500D Cordura® solution dyed nylon
- Sand- and water-resistant, self repairing zippers with silent 550 pull cords
- Separate lid compartment to secure the R.E.E.L. Non-Rigid Litter
- 2 internal lid pouches to secure additional medical products
- 2 external side pouches on either side of the kit for additional medical products
- 2 external straps with dual adjust MILSPEC IR reduced buckles to secure the kit in a vehicle or boat

ECRK

ITEM #	NSN	
80-0889	----	R _X ONLY

* Explore other trauma management kits at www.NARescue.com

DIMENSIONS


- Case closed: H 17 in. x W 16.5 in. x D 10 in.
- Weight: 18 lb



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

EXPEDITIONARY CASUALTY RESPONSE KIT (ECRK)

KIT INCLUDES:



KIT CONTENTS

- 6 x Combat Application Tourniquet® (C-A-T®) - Black
- 1 x Blizzard Blanket
- 4 x NAR® Nasopharyngeal Airway with Lube - 28F(117 mm)
- 4 x ARS® Needle Decompression Device, 14g x 3.25 in.
- 4 x Combat Gauze, Z-Fold Hemostatic Dressing
- 6 x Emergency Trauma Dressing (ETD®) 6 in.
- 10 x NAR® Compressed Gauze, 4.5 in. x 4.1 yd
- 4 x USMC HyFin® Chest Seal Combo Pk (1 HyFin® Vent & 1 HyFin® Chest Seal)
- 6pr Bear Claw™ Gloves, Large Tan
- 6 x Emergency Trauma Dressing (ETD®) ABD 8 x 10 in.
- 1 x Minor BurnTec® Kit
- 1 x Trauma Shears, Large, 7.25 in.
- 6 x Combat Casualty Response Cards
- 2 x Surgical Tape 1 in.
- 1 x R.E.E.L. Non-Rigid Litter
- 10 x Antiseptic Wipes
- 1 x Combat Casualty Reference Card
- 2 x SAM® Splint II
- 2 x Large Permanent Marker, Black
- 3 x Elastic Bandage 6 in.
- 2 x Triangular Bandage, 45 in. x 45 in. x 63 in.
- 2 x Mini Duct Tape, 2 in. x 100 in.
- 1 x Model 7 Rescue Hook
- 6 x Polycarbonate Eye Shield
- 1 x Expeditionary Casualty Response Bag



Learn more about this and other Products with a Mission® at www.NARescue.com



ECRK

ITEM #	NSN	
80-0889	----	Rx ONLY

* Explore other trauma management kits at www.NARescue.com