

RUF KIT

Slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with additional trauma supplies

The RUF Kit from North American Rescue is a slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with additional trauma supplies. Designed with two tear away internal panels that hold standard IFAK contents and two tear away clear accessory pouches for additional supplies. The internal kit also has a shock cord panel designed to hold a compact, soft litter. The kit can be worn with a single bandolier style sling strap as well as comes with a waist strap that can be removed based on the user's needs.

The RUF Bag is also available as a bag only (No gear or PVC Med Patch)



Learn more about this and other Products with a Mission® at www.NARescue.com



FEATURES/SPECIFICATIONS

- Slim profile trauma aid pack with two sets of IFAK components plus extra trauma supplies
- Can be worn with single bandolier style sling strap and with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Internal shock cord area to hold compact soft litter or additional equipment
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

RUF KIT

ITEM #	OPTION
80-0974	Basic
80-0975	w/ Combat Gauze LE
80-0976	Advanced
80-0967W	Bag only

DIMENSIONS

- Dimensions: H 18 in. x W 8 in. x D 2 in.
- Kit Weight: 6.3 lb.
- Bag Weight: 2 lb.




NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

RUF KIT

Every RUF KIT INCLUDES THESE CORE ITEMS:



+plus, by kit type:



KIT CONTENTS

- 2 x C-A-T® Tourniquet - Black
- 2 x NAR® S-Rolled Gauze 4.5 in. x 4.1 yd.

+plus, by kit type:

- 2 x NAR® S-Rolled Gauze - Or -
- 2 x Combat Gauze® LE - Or -
- 2 x Combat Gauze® LE
- 2 x NAR® Nasopharyngeal Airway with Lube, 28F
- 2 x ARS® Needle - 14g

- 2 x ETD® Dressing - 6 in.
- 1 x ETD® Abdominal / Stump
- 2 x Hyfin Vent Compact Twin Pack
- 2 x Emergency Survival Wrap - Orange
- 1 x NAR® ELITE™ - LITE Responder Litter
- 4pr Nitrile Gloves - Large
- 1 x Responder Shears - Large
- 2 x Permanent Marker
- 1 x Black/Gray Patch - 2 in. x 2 in.



Can be worn with single bandolier style sling strap as well as with a detachable waist strap



Learn more about this and other Products with a Mission® at www.NARescue.com

